

Biologics Manufacturing & *GMP Quality Control*

Twenty-five industrial training programmes, one for each distinct type of NGS data analysis performed in biologics, cell and gene therapy manufacturing under quality systems — from cell bank characterisation and viral safety through vector and product release testing to method validation and data integrity. Every programme answers the questions a laboratory asks before it commits budget — why the analysis is done, when it is done, what it detects, what it does not detect, what makes it hard, who does this job and what a trainee can do afterwards — and is followed by its complete ordered workflow, step by step.

25

ANALYSIS TYPES

481

STEPS

2328

HOURS

397

DAYS

Biologics Manufacturing & GMP Quality Control — *Analysis Types*

Twenty-five distinct types of NGS data analysis performed in regulated biologics manufacturing, each an independent analysis with its own input, its own pipeline and its own report set. Each is listed under its formal title with the industry name and the shorthand a laboratory uses in practice. Every analysis type occupies a two-page entry: the first page answers why, when, what it detects, what it does not detect, who does the job and what follows; the second sets out the complete ordered workflow.

CODE	ANALYSIS TYPE · INDUSTRY NAME & SHORT NAME	STEPS	HOURS	DAYS	FEE (INR)
SBQ-01	Master Cell Bank Genomic Characterisation Master cell bank genomic characterisation · MCB characterisation	22	106	18	2,45,000
SBQ-02	Working Cell Bank & End-of-Production Cell Analysis WCB and EOPC genomic testing · WCB · EOPC analysis	19	92	16	2,10,000
SBQ-03	Cell Line Genetic Stability Analysis Producer cell line genetic stability assessment · Genetic stability analysis	20	96	16	2,20,000
SBQ-04	Cell Substrate Identity & Purity Analysis Cell substrate identity and cross-contamination testing · Substrate identity testing	16	76	13	1,75,000
SBQ-05	Adventitious Virus Detection by Sequencing Sequencing-based adventitious virus testing · Viral safety sequencing	23	112	19	2,55,000
SBQ-06	Bacterial & Mycoplasma Detection by Sequencing Sequencing-based microbial and mycoplasma testing · Microbial detection sequencing	18	86	15	2,00,000
SBQ-07	Endogenous Retroviral Element Analysis Endogenous retroviral element characterisation · ERV analysis	18	88	15	2,05,000
SBQ-08	Raw Material & Feedstock Screening Analysis Raw material adventitious agent screening · Raw material screening	17	82	14	1,90,000
SBQ-09	Transgene Integrity & Copy Number Analysis Transgene integration and copy number characterisation · Transgene characterisation	19	92	16	2,10,000
SBQ-10	Expression Construct Verification Analysis Expression construct sequence verification · Construct verification	17	82	14	1,90,000
SBQ-11	Producer Cell Line Clone Selection Analysis Clone selection genomic screening · Clone selection analysis	20	96	16	2,20,000
SBQ-12	Host Cell DNA Residual Analysis Residual host cell DNA characterisation · Residual DNA analysis	18	86	15	2,00,000
SBQ-13	Viral Vector Genome Integrity Analysis Vector genome integrity and impurity characterisation · Vector integrity QC	22	106	18	2,45,000

Continued overleaf. Programme codes read SBQ-01 and protocol references PRT/SBQ/01.

Analysis Types — 14 to 25

CODE	ANALYSIS TYPE · INDUSTRY NAME & SHORT NAME	STEPS	HOURS	DAYS	FEE (INR)
SBQ-14	Capsid & Packaging Composition Analysis Vector packaging composition characterisation · Packaging composition analysis	20	98	17	2,25,000
SBQ-15	Plasmid Identity & Quality Analysis Plasmid sequence identity and quality testing · Plasmid QC	17	82	14	1,90,000
SBQ-16	mRNA Sequence & Impurity Analysis mRNA product sequence and impurity characterisation · mRNA QC	20	96	16	2,20,000
SBQ-17	Cell Therapy Product Release Analysis Cell therapy genomic release testing · Cell therapy release testing	22	108	18	2,45,000
SBQ-18	Microbial Identification for Contamination Investigation Contamination isolate identification and source tracing · Contamination investigation	18	86	15	2,00,000
SBQ-19	Environmental & Utility Monitoring Metagenomics Facility environmental monitoring by sequencing · Environmental monitoring metagenomics	17	82	14	1,90,000
SBQ-20	Process Comparability Genomic Analysis Manufacturing comparability genomic assessment · Comparability analysis	21	102	17	2,35,000
SBQ-21	Stability Study Genomic Analysis Product stability genomic assessment · Stability genomic analysis	16	78	13	1,80,000
SBQ-22	Analytical Method Validation for Sequencing Assays Sequencing assay validation for GMP use · Method validation	24	118	20	2,70,000
SBQ-23	Data Integrity & Computerised System Compliance Analysis Data integrity assessment for computational systems · Data integrity compliance	20	98	17	2,25,000
SBQ-24	Deviation & Out-of-Specification Investigation Analysis Genomic investigation of manufacturing deviations · Deviation investigation	19	92	16	2,10,000
SBQ-25	Rapid Sterility Assurance by Sequencing Methods Rapid microbiological method implementation · Rapid sterility methods	18	88	15	2,05,000
Complete Sector Catalogue 02 — all thirty-nine analysis types		481	2328	397	53,60,000

Each fee is derived independently from the step count, analytical complexity and reporting burden of that specific analysis type. There is no banded pricing. Days are calculated at six contact hours per working day and rounded up. Fees are per participant and exclusive of GST and the one-time registration charge. In-house delivery at the client laboratory and institutional rates are available on application.

WHY IT IS DONE

The master cell bank is the permanent origin of every batch a product will ever produce, so it is characterised once, exhaustively, and never again from the same material. Anything missed at this point propagates through the entire commercial life of the product.

WHAT IT ANSWERS

What is the complete genomic identity, integrity and purity status of this master cell bank, and does it support use as the origin of an approved product?

WHEN IT IS DONE

At master cell bank establishment before any further banking, and when a legacy bank requires re-characterisation to modern expectations.

WHAT TRIGGERS IT

Establishment of a new master cell bank, a regulatory requirement to re-characterise a legacy bank, or a technology transfer requiring full characterisation.

HOW IT IS DONE

Identity and species are established, the host genome and integrated construct are sequenced and characterised, and integration sites and copy number are determined. Adventitious agent screening runs against comprehensive references with controls processed identically, endogenous retroviral elements are characterised, and the complete profile is archived as the permanent baseline.

WHAT IT PRODUCES

Identity and species confirmation, host genome characterisation, construct sequence verification, integration site and copy number determination, adventitious agent screening results with controls, endogenous element characterisation, and the archived baseline profile.

WHAT IT DETECTS

Cell line identity and species, genomic sequence of the host and any integrated construct, transgene integration sites and copy number, adventitious agent sequence, endogenous retroviral element status, and the baseline genomic state against which all later banks are compared.

WHAT IT DOES NOT DETECT

It does not establish productivity, product quality attributes or viability, which are separate release attributes. Nucleic acid detection of an agent does not establish infectivity, and genomic characterisation cannot anticipate every process-related change that will occur downstream.

WHAT MAKES IT HARD

This characterisation is done once and referenced for decades, so anything not captured cannot be recovered without re-establishing the bank. Distinguishing genuine agent detection from ubiquitous reagent background requires controls designed in from the start rather than added after a positive appears.

WHO DOES THIS JOB

Bioinformatician in biologics quality or cell bank characterisation, senior level, working with quality assurance and regulatory affairs.

WHAT YOU CAN DO AFTERWARDS

Characterise a master cell bank to regulatory expectations, establish a permanent genomic baseline, and produce documentation that supports a marketing application.

WHO HIRES FOR IT

Biologics manufacturers, cell and gene therapy companies, contract development and manufacturing organisations, and cell bank testing providers.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
22	106 hours · 18 days	Prior genome analysis experience, GMP awareness	2,45,000

COMPLETE WORKFLOW — 22 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Characterisation protocol authorisation	The characterisation protocol and its regulatory basis approved before work begins	Protocol approval	<i>Authorised protocol</i>
02	Bank provenance documentation	Cell line origin, derivation history and all prior handling documented	Provenance records	<i>Provenance record</i>
03	Sample chain of custody establishment	Custody from bank vial to analytical sample established and recorded	Custody procedure	<i>Custody record</i>
04	Species identification	Species confirmed independently of any claimed identity	Species analysis	<i>Species confirmation</i>
05	Identity profile generation	An identity profile generated and archived as the permanent reference	Genotyping	<i>Identity profile</i>
06	Data receipt and quality assessment	Sequence integrity and quality verified against acceptance criteria	md5sum, FastQC	<i>QC acceptance record</i>
07	Host genome alignment and processing	Reads aligned to the host reference and processed for analysis	bwa-mem2, Picard	<i>Analysis-ready BAM</i>
08	Genome coverage assessment	Coverage and uniformity assessed against characterisation requirements	mosdepth	<i>Coverage record</i>
09	Host genome variant characterisation	Genomic variation relative to the reference characterised as the baseline	Variant calling	<i>Host variant profile</i>
10	Karyotypic and copy number assessment	Chromosomal complement and copy number state characterised	Copy number analysis	<i>Karyotypic profile</i>
11	Construct sequence verification	Any integrated construct verified against its design sequence	Sequence comparison	<i>Construct verification</i>
12	Integration site mapping	Integration sites mapped to host coordinates with junctions recovered	Junction analysis	<i>Integration sites</i>
13	Transgene copy number quantification	Copy number quantified against host genome background	Copy number quantification	<i>Copy number result</i>
14	Integration structure reconstruction	Integration architecture reconstructed including concatemeric arrangements	Structure reconstruction	<i>Integration structure</i>
15	Host depletion for agent screening	Host sequence depleted with efficiency quantified before agent classification	Host depletion	<i>Non-host reads</i>
16	Adventitious agent classification	Non-host reads classified against comprehensive reference databases	Classification	<i>Agent classification</i>
17	Control comparison	Detections assessed against reagent and negative controls processed identically	Control analysis	<i>Control comparison</i>
18	Spike-in sensitivity determination	Achieved detection sensitivity established by spike-in across agent classes	Spike analysis	<i>Sensitivity determination</i>
19	Endogenous retroviral element characterisation	Endogenous elements characterised and distinguished from exogenous sequence	Element analysis	<i>Endogenous profile</i>
20	Unassigned sequence characterisation	Sequence matching no reference assembled and characterised	Assembly, characterisation	<i>Unassigned characterisation</i>
21	Baseline profile archiving	The complete characterisation archived as the permanent reference baseline	Archive procedure	<i>Archived baseline</i>
22	Characterisation report assembly	All results assembled into a regulatory characterisation dossier	Report assembly	<i>Characterisation dossier</i>

STEPS IN WORKFLOW

22

TRAINING DURATION

106 h · 18 d

PROGRAMME CODE

SBQ-01

TRAINING FEE (INR)

2,45,000

Working Cell Bank & End-of-Production Cell Analysis

(WCB and EOPC genomic testing · WCB · EOPC analysis)

WHY IT IS DONE

Cells change through expansion, and the population at the end of production has been through many more doublings than the bank it came from. Testing at both points brackets the manufacturing window and demonstrates that the process does not drive the cells outside their characterised state.

WHAT IT ANSWERS

Do the working cell bank and the end-of-production cells remain consistent with the master bank, and has anything changed across the manufacturing window?

WHEN IT IS DONE

At working cell bank establishment, at the limit of in vitro cell age during process characterisation, and at end of production for each qualified process.

WHAT TRIGGERS IT

Working cell bank establishment, process characterisation studies, or a limit of in vitro cell age qualification requirement.

HOW IT IS DONE

Both banks and end-of-production cells are analysed against the archived master bank baseline, identity concordance is confirmed, and genomic differences are identified and assessed. Transgene copy number and integrity are quantified at each point, adventitious agent screening is repeated, and changes are characterised across the manufacturing window.

WHAT IT PRODUCES

Identity concordance against the master bank, genomic change assessment across expansion, transgene copy number and integrity at each point, adventitious agent screening results, and a manufacturing window stability statement.

WHAT IT DETECTS

Identity concordance with the master bank, genomic changes accumulated through expansion, transgene copy number and integrity at the production limit, adventitious agent status at both points, and any population shift relative to the bank.

WHAT IT DOES NOT DETECT

It cannot establish that changes observed are inconsequential for product quality, which requires product characterisation. Low-frequency subpopulations below detection may still expand under different process conditions.

WHAT MAKES IT HARD

Deciding which observed changes matter is the difficulty, since some genomic drift is inevitable and not all of it is consequential. The comparison depends entirely on the master bank baseline having been characterised well enough to compare against.

WHO DOES THIS JOB

Bioinformatician in biologics quality control, mid to senior level, working with process development.

WHAT YOU CAN DO AFTERWARDS

Compare banks and production cells against a characterised baseline, assess genomic change across the manufacturing window, and support limit of in vitro cell age qualification.

WHO HIRES FOR IT

Biologics manufacturers, contract manufacturing organisations, cell and gene therapy companies, and quality control testing laboratories.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
19	92 hours · 16 days	Prior cell bank characterisation experience	2,10,000

Working Cell Bank & End-of-Production Cell Analysis

(WCB and EOPC genomic testing · WCB · EOPC analysis)

COMPLETE WORKFLOW — 19 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Testing protocol confirmation	The applicable testing protocol and acceptance criteria confirmed	Protocol review	Testing protocol
02	Sample provenance recording	Bank or production stage, passage and custody recorded for each sample	Sample records	Provenance record
03	Master bank baseline retrieval	The archived master bank baseline retrieved for comparison	Archive retrieval	Retrieved baseline
04	Data receipt and quality assessment	Sequence integrity and quality verified per sample	md5sum, FastQC	QC acceptance record
05	Alignment and processing	Reads aligned and processed identically to the baseline characterisation	bwa-mem2, Picard	Analysis-ready BAM
06	Identity profile generation	Identity profile generated for concordance against the master bank	Genotyping	Identity profile
07	Identity concordance assessment	Concordance with the master bank baseline assessed against criteria	Concordance analysis	Concordance determination
08	Coverage assessment	Coverage assessed against the requirements of each comparison	mosdepth	Coverage record
09	Genomic difference identification	Differences from the baseline identified across variant and copy number	Comparative analysis	Difference set
10	Difference significance assessment	Identified differences assessed for likely consequence	Assessment framework	Significance assessment
11	Transgene copy number quantification	Copy number quantified and compared against the baseline	Copy number quantification	Copy number comparison
12	Transgene integrity verification	Construct sequence integrity verified at the production stage	Sequence comparison	Integrity verification
13	Integration site verification	Integration sites confirmed unchanged from the baseline	Junction analysis	Site verification
14	Host depletion and agent screening	Host depleted and adventitious agent screening performed	Depletion, classification	Agent screening results
15	Control comparison	Detections assessed against controls processed identically	Control analysis	Control comparison
16	Population shift assessment	Clonal composition compared against the baseline for population shift	Clonal analysis	Population assessment
17	Manufacturing window assessment	Findings assessed across the full window from bank to production limit	Window analysis	Window assessment
18	Acceptance determination	Results assessed against protocol acceptance criteria	Acceptance criteria	Acceptance determination
19	Testing report assembly	Concordance, differences and agent screening assembled into the testing report	Report assembly	Testing report

STEPS IN WORKFLOW	TRAINING DURATION	PROGRAMME CODE	TRAINING FEE (INR)
19	92 h · 16 d	SBQ-02	2,10,000

Cell Line Genetic Stability Analysis

(Producer cell line genetic stability assessment · Genetic stability analysis)

WHY IT IS DONE

A production cell line that loses transgene copies, rearranges its integration site or accumulates chromosomal change over expansion will eventually produce a different product from the one that was approved. Stability data define how long a line can be used.

WHAT IT ANSWERS

Is this cell line genetically stable across the intended manufacturing window, and where does instability begin?

WHEN IT IS DONE

In process development to establish the limit of in vitro cell age, during process qualification, and when investigating a productivity or product quality drift.

WHAT TRIGGERS IT

Establishment of a production cell line, a limit of in vitro cell age study, or an unexplained drift in productivity or product quality.

HOW IT IS DONE

Cells are sampled at defined passage intervals across and beyond the intended window, transgene copy number is quantified at each, and integration sites are re-characterised to detect rearrangement. Chromosomal change is assessed, clonal composition is tracked to detect population shift, and the passage at which change appears is identified.

WHAT IT PRODUCES

Copy number trajectory across passage, integration site integrity at each interval, chromosomal change assessment, clonal composition tracking, the detected onset of instability, and a limit of in vitro cell age recommendation.

WHAT IT DETECTS

Transgene copy number change across passage, integration site rearrangement or loss, chromosomal-level change including aneuploidy, clonal population shift within the culture, and the passage at which any change becomes detectable.

WHAT IT DOES NOT DETECT

It does not link genomic change to product quality consequence, which requires parallel product characterisation. Changes affecting only a minority subpopulation may be below detection until that population expands.

WHAT MAKES IT HARD

Population-level measurement conceals subclonal change until it dominates, so tracking clonal composition rather than only average copy number is what detects instability early. Sampling intervals must extend beyond the intended window or the study cannot demonstrate a margin.

WHO DOES THIS JOB

Bioinformatician in process development or quality, mid to senior level.

WHAT YOU CAN DO AFTERWARDS

Design and analyse genetic stability studies, track clonal composition rather than averages, and derive defensible limits of in vitro cell age.

WHO HIRES FOR IT

Biologics manufacturers, contract development organisations, cell line development companies, and process development functions.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
20	96 hours · 16 days	Prior copy number and structural variant experience	2,20,000

COMPLETE WORKFLOW — 20 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Study design authorisation	Passage intervals, endpoints and the intended window approved before study start	Protocol approval	Authorised design
02	Passage interval definition	Sampling intervals defined to extend beyond the intended manufacturing window	Design specification	Sampling schedule
03	Culture and sampling record capture	Culture conditions and sampling events recorded for every timepoint	Culture records	Culture record
04	Baseline characterisation	The starting population characterised as the comparison baseline	Baseline analysis	Baseline profile
05	Data receipt and quality assessment	Sequence integrity and quality verified per timepoint	md5sum, FastQC	QC acceptance record
06	Alignment and processing	All timepoints aligned and processed identically	bwa-mem2, Picard	Analysis-ready BAMs
07	Coverage consistency assessment	Coverage confirmed consistent across timepoints for valid comparison	mosdepth	Coverage consistency
08	Transgene copy number quantification	Copy number quantified at every passage interval	Copy number quantification	Copy number trajectory
09	Copy number trend modelling	Copy number modelled as a trend across passage rather than pairwise	Trend modelling	Copy number model
10	Integration junction verification	Integration junctions re-verified at each interval for rearrangement	Junction analysis	Junction verification
11	Integration structure comparison	Integration architecture compared across intervals	Structure comparison	Structure comparison
12	Transgene sequence integrity assessment	Construct sequence integrity assessed at each interval	Sequence comparison	Integrity assessment
13	Chromosomal change assessment	Chromosomal complement and copy number change assessed across passage	Copy number analysis	Chromosomal assessment
14	Clonal composition analysis	Clonal composition characterised at each interval	Clonal analysis	Clonal composition
15	Population shift detection	Shifts in dominant clonal populations detected across the study	Shift analysis	Population shift
16	Subclonal emergence assessment	Emerging subpopulations detected before they dominate	Subclonal analysis	Emergence assessment
17	Instability onset determination	The passage at which any change becomes detectable determined	Onset analysis	Onset determination
18	Margin assessment	Margin between instability onset and the intended window assessed	Margin analysis	Margin assessment
19	Limit of in vitro cell age recommendation	A defensible limit recommended from the observed trajectories	Derivation framework	Limit recommendation
20	Stability report assembly	Trajectories, onset and limit recommendation assembled	Report assembly	Stability report

STEPS IN WORKFLOW	TRAINING DURATION	PROGRAMME CODE	TRAINING FEE (INR)
20	96 h · 16 d	SBQ-03	2,20,000

Cell Substrate Identity & Purity Analysis

(Cell substrate identity and cross-contamination testing · Substrate identity testing)

WHY IT IS DONE

A manufacturing facility running multiple cell lines carries a permanent risk of cross-contamination between them, and a contaminated substrate produces a product that is not what the marketing authorisation describes. Identity and purity testing is the control against that.

WHAT IT ANSWERS

Is this cell substrate the declared line, is it free of other cell lines, and does it match the reference profile held for it?

WHEN IT IS DONE

On receipt of any substrate, at bank establishment, at defined intervals during campaigns, and in investigating a suspected cross-contamination event.

WHAT TRIGGERS IT

Substrate receipt, bank establishment, a scheduled interval, or a contamination investigation.

HOW IT IS DONE

An identity profile is generated and compared against the archived reference for that substrate, species is confirmed independently, and mixed allele signal is analysed for evidence of a second line. Contamination proportion is estimated where present, and the facility's other lines are compared as candidate sources.

WHAT IT PRODUCES

Identity profile with reference concordance, species confirmation, cross-contamination assessment with proportion, candidate source identification from facility lines, and a purity determination.

WHAT IT DETECTS

Identity profile against the reference, species confirmation, cross-contamination from other cell lines with proportion estimated, and inter-species contamination.

WHAT IT DOES NOT DETECT

It cannot detect contamination by a line with an identical profile, and it does not assess microbiological contamination, which requires separate methods. Very low-level cross-contamination falls below the discrimination limit of any marker panel.

WHAT MAKES IT HARD

Identifying the source line is what converts a detection into a corrective action, and that requires reference profiles for every line the facility handles. Low-level contamination sits near the discrimination limit and must be reported with its detection bound.

WHO DOES THIS JOB

Bioinformatician or quality control analyst, entry to mid level, within a quality function.

WHAT YOU CAN DO AFTERWARDS

Perform substrate identity and purity testing, estimate contamination proportion, and identify candidate source lines within a facility.

WHO HIRES FOR IT

Biologics manufacturers, contract manufacturing organisations, cell banking services, and quality control laboratories.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
16	76 hours · 13 days	Prior genotype analysis experience	1,75,000

Cell Substrate Identity & Purity Analysis

(Cell substrate identity and cross-contamination testing · Substrate identity testing)

COMPLETE WORKFLOW — 16 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Testing scope confirmation	The substrate, testing occasion and applicable criteria confirmed	Protocol review	<i>Testing scope</i>
02	Reference profile retrieval	The archived reference profile for the substrate retrieved	Archive retrieval	<i>Reference profile</i>
03	Facility line inventory retrieval	Profiles of all other lines handled in the facility retrieved for comparison	Archive retrieval	<i>Facility profile set</i>
04	Sample custody verification	Custody from culture to analytical sample verified	Custody procedure	<i>Custody record</i>
05	Data receipt and quality assessment	Data quality verified against acceptance criteria	Quality assessment	<i>QC acceptance record</i>
06	Identity profile generation	The identity profile generated across the marker panel	Genotyping	<i>Identity profile</i>
07	Profile completeness assessment	Profile completeness and call quality assessed	Quality analysis	<i>Profile quality</i>
08	Reference concordance assessment	Profile compared against the archived reference with thresholds applied	Concordance analysis	<i>Concordance determination</i>
09	Species confirmation	Species confirmed independently of the identity match	Species analysis	<i>Species confirmation</i>
10	Mixed allele analysis	Additional alleles indicating a second cell line identified	Mixture analysis	<i>Mixture findings</i>
11	Contamination proportion estimation	Where mixed, the proportion of the contaminant estimated with bounds	Proportion estimation	<i>Contamination estimate</i>
12	Facility line comparison	Mixed signal compared against facility line profiles to identify the source	Comparative analysis	<i>Source identification</i>
13	Detection limit statement	The contamination detection limit of the marker panel stated	Limit calculation	<i>Detection limit</i>
14	Inter-species contamination assessment	Contamination from other species assessed separately	Species analysis	<i>Inter-species assessment</i>
15	Purity determination	A purity determination made against the applicable criteria	Determination criteria	<i>Purity determination</i>
16	Testing report assembly	Identity, purity, source and detection limit assembled	Report assembly	<i>Identity and purity report</i>

STEPS IN WORKFLOW

16

TRAINING DURATION

76 h · 13 d

PROGRAMME CODE

SBQ-04

TRAINING FEE (INR)

1,75,000

Adventitious Virus Detection by Sequencing

(Sequencing-based adventitious virus testing · Viral safety sequencing)

WHY IT IS DONE

Conventional viral safety testing detects only what its assays are designed for, and historical contamination events involved viruses those assays could not see. Sequencing detects any viral sequence present, which is why regulators now accept it as a broad complement to targeted testing.

WHAT IT ANSWERS

Does this material contain any viral sequence, of what identity, at what abundance, and does the evidence support release?

WHEN IT IS DONE

In cell bank and end-of-production testing, in bulk harvest testing, in raw material qualification, and in investigating a suspected contamination event.

WHAT TRIGGERS IT

Cell bank or harvest testing, raw material qualification, a suspect observation in cell culture, or a contamination investigation.

HOW IT IS DONE

Host sequence is depleted with efficiency quantified because it sets the effective sensitivity, non-host reads are classified against comprehensive viral references, and every detection is assessed against reagent and negative controls processed identically. Spike-in controls establish achieved sensitivity, detections are confirmed by targeted analysis, and unassigned sequence is assembled and characterised.

WHAT IT PRODUCES

Host depletion efficiency, viral detections with identity and abundance, control comparison establishing background, spike-in derived sensitivity per agent class, targeted confirmation of detections, unassigned sequence characterisation, and a viral safety assessment.

WHAT IT DETECTS

Viral sequence from any family present above the detection threshold, species and strain identification where references exist, viral abundance relative to host, integrated and episomal viral sequence, and novel viral sequence flagged as unassigned.

WHAT IT DOES NOT DETECT

It detects nucleic acid rather than infectious virus, so significance requires follow-up. Viruses absent from reference databases can only be flagged as unassigned sequence, and detection sensitivity depends entirely on abundance relative to the overwhelming host background.

WHAT MAKES IT HARD

Sensitivity must be demonstrated rather than claimed, and it varies by an order of magnitude between agents depending on genome size and host background. Reagent-borne viral sequence is ubiquitous, so controls processed in parallel are what separate a finding from laboratory background.

WHO DOES THIS JOB

Bioinformatician in viral safety or biosafety testing, senior level, working with virology and quality assurance.

WHAT YOU CAN DO AFTERWARDS

Run sequencing-based viral safety testing with demonstrated sensitivity, distinguish findings from reagent background, and produce regulatory-grade safety assessments.

WHO HIRES FOR IT

Biologics manufacturers, viral safety testing laboratories, cell and gene therapy companies, and contract testing organisations.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
23	112 hours · 19 days	Prior metagenomic analysis experience	2,55,000

Adventitious Virus Detection by Sequencing

(Sequencing-based adventitious virus testing · Viral safety sequencing)

COMPLETE WORKFLOW — 23 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Testing protocol authorisation	Scope, agent classes and acceptance criteria approved before testing	Protocol approval	Authorised protocol
02	Reference database assembly	Comprehensive viral reference databases assembled with versions recorded	Database assembly	Reference databases
03	Control design	Negative, reagent and spike-in controls designed into the run	Control design	Control plan
04	Spike-in panel preparation	Spike organisms spanning genome types prepared at defined levels	Spike preparation	Spike controls
05	Sample and control extraction	Samples and controls extracted in parallel under identical conditions	Extraction	Prepared material
06	Nucleic acid type coverage confirmation	Both DNA and RNA workflows confirmed where the scope requires both	Workflow verification	Coverage confirmation
07	Data receipt and quality assessment	Sequence integrity and run quality verified per sample and control	md5sum, FastQC	QC acceptance record
08	Host sequence depletion	Host reads removed with the method recorded	Host depletion	Non-host reads
09	Depletion efficiency quantification	Depletion efficiency quantified as it determines effective sensitivity	Efficiency analysis	Depletion metrics
10	Viral classification	Non-host reads classified against the viral reference databases	Classification	Viral classification
11	Read distribution assessment	Read distribution across each detected genome assessed for authenticity	Coverage analysis	Distribution assessment
12	Reagent background subtraction	Taxa present in reagent controls treated as background	Background analysis	Background-corrected results
13	Negative control comparison	Sample detections compared against parallel negative controls	Control comparison	Control comparison record
14	Spike-in recovery quantification	Spike organism recovery quantified per genome type	Recovery analysis	Recovery metrics
15	Per-class sensitivity determination	Detection sensitivity determined per viral genome class	Sensitivity calculation	Per-class sensitivity
16	Detection threshold application	Detections assessed against the reporting threshold	Threshold criteria	Threshold decisions
17	Targeted confirmation	Detections confirmed by targeted analysis or genome assembly	Confirmation analysis	Confirmation results
18	Integrated viral sequence assessment	Viral sequence integrated into the host genome identified separately	Integration analysis	Integration findings
19	Endogenous element discrimination	Endogenous retroviral sequence distinguished from exogenous virus	Discrimination analysis	Endogenous discrimination
20	Unassigned sequence assembly	Unclassified sequence assembled and characterised	Assembly, characterisation	Unassigned characterisation
21	Viability limitation documentation	The distinction between nucleic acid and infectious virus documented	Limitation template	Viability caveat
22	Risk assessment	Detections assessed for significance against the material's intended use	Risk assessment	Risk determination
23	Viral safety report assembly	Detections, sensitivity, controls and risk assembled into a safety report	Report assembly	Viral safety report

STEPS IN WORKFLOW

23

TRAINING DURATION

112 h · 19 d

PROGRAMME CODE

SBQ-05

TRAINING FEE (INR)

2,55,000

Bacterial & Mycoplasma Detection by Sequencing

(Sequencing-based microbial and mycoplasma testing · Microbial detection sequencing)

WHY IT IS DONE

Mycoplasma contamination is common, invisible in culture and alters cell behaviour profoundly. Sequencing detects it and other bacterial contaminants without culture, in days rather than the weeks compendial methods require, which changes what a manufacturer can do about it.

WHAT IT ANSWERS

Does this material contain bacterial or mycoplasmal sequence, of what identity, and does it represent contamination or reagent background?

WHEN IT IS DONE

In routine cell culture monitoring, in bank and harvest testing, in raw material qualification, and in investigating a suspected contamination.

WHAT TRIGGERS IT

Routine monitoring, bank or harvest testing, an anomalous culture observation, or a contamination investigation.

HOW IT IS DONE

Host sequence is depleted, non-host reads are classified against curated bacterial and mycoplasmal references, and detections are assessed against reagent and negative controls run in parallel. Species assignment is confirmed by marker gene analysis or assembly, abundance is quantified, and spike-in controls establish sensitivity.

WHAT IT PRODUCES

Classification results with species assignment, abundance quantification, control comparison establishing reagent background, spike-in derived sensitivity, confirmation analysis for detections, and a contamination determination.

WHAT IT DETECTS

Bacterial and mycoplasmal sequence with species-level identification where references permit, relative abundance, and the distinction between sample contamination and reagent-derived background.

WHAT IT DOES NOT DETECT

It detects nucleic acid rather than viable organism, so a positive does not establish live contamination and a negative does not establish sterility. Species absent from references are identifiable only to higher taxonomic rank.

WHAT MAKES IT HARD

Reagent and environmental bacterial sequence is present in every preparation, so distinguishing genuine contamination from background is the entire analytical problem and requires controls rather than thresholds alone. Low-biomass samples are dominated by background.

WHO DOES THIS JOB

Bioinformatician in quality control or biosafety, mid level.

WHAT YOU CAN DO AFTERWARDS

Detect and identify microbial contamination by sequencing, distinguish it from reagent background, and demonstrate sensitivity by spike-in.

WHO HIRES FOR IT

Biologics manufacturers, quality control laboratories, contract testing organisations, and cell therapy manufacturers.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
18	86 hours · 15 days	Prior metagenomic analysis experience	2,00,000

Bacterial & Mycoplasma Detection by Sequencing

(Sequencing-based microbial and mycoplasma testing · Microbial detection sequencing)

COMPLETE WORKFLOW — 18 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Testing scope confirmation	Material type, agent scope and acceptance criteria confirmed	Protocol review	Testing scope
02	Reference database assembly	Bacterial and mycoplasmal references assembled with versions recorded	Database assembly	Reference databases
03	Control design and preparation	Negative, reagent and spike controls designed and prepared	Control preparation	Control set
04	Sample and control extraction	Samples and controls extracted in parallel under identical conditions	Extraction	Prepared material
05	Data receipt and quality assessment	Sequence integrity and quality verified per sample and control	md5sum, FastQC	QC acceptance record
06	Host sequence depletion	Host reads removed with efficiency quantified	Host depletion	Non-host reads
07	Taxonomic classification	Non-host reads classified against bacterial and mycoplasmal references	Classification	Classification results
08	Marker gene analysis	Detections confirmed and refined by marker gene analysis	Marker analysis	Marker confirmation
09	Reagent background subtraction	Taxa present in reagent controls treated as background	Background analysis	Corrected results
10	Negative control comparison	Detections compared against parallel negative controls	Control comparison	Control comparison record
11	Spike-in recovery assessment	Spike organism recovery quantified to establish sensitivity	Recovery analysis	Sensitivity determination
12	Abundance quantification	Detected organism abundance quantified relative to host and total	Quantification	Abundance results
13	Species assignment confidence	Assignment confidence assessed with higher-rank assignment where uncertain	Confidence assessment	Assignment confidence
14	Threshold application	Detections assessed against the reporting abundance threshold	Threshold criteria	Threshold decisions
15	Genome assembly for confirmation	Detected organisms assembled where abundance permits confirmation	Assembly	Assembly confirmation
16	Viability limitation documentation	The nucleic acid versus viable organism distinction documented	Limitation template	Viability caveat
17	Contamination determination	A contamination determination made against acceptance criteria	Determination criteria	Contamination determination
18	Testing report assembly	Detections, sensitivity, controls and determination assembled	Report assembly	Microbial testing report

STEPS IN WORKFLOW

18

TRAINING DURATION

86 h · 15 d

PROGRAMME CODE

SBQ-06

TRAINING FEE (INR)

2,00,000

Endogenous Retroviral Element Analysis

(Endogenous retroviral element characterisation · ERV analysis)

WHY IT IS DONE

Rodent-derived production cell lines carry endogenous retroviral sequences that are expressed and can produce particles. Characterising them and demonstrating that downstream processing clears them is a defined requirement of viral safety for these products.

WHAT IT ANSWERS

What endogenous retroviral elements does this cell substrate carry, are they expressed, and what is their contribution to the particle burden the process must clear?

WHEN IT IS DONE

In master cell bank characterisation, in process viral clearance justification, and in supporting the viral safety section of a regulatory submission.

WHAT TRIGGERS IT

Master cell bank characterisation, a viral safety assessment, or a regulatory question on endogenous element burden.

HOW IT IS DONE

Host genome sequence is analysed for retroviral element copies with positions mapped, element structure is assessed for intact reading frames, and expression is quantified from transcriptomic data. Elements are distinguished from any exogenous retroviral sequence, and burden is quantified for viral clearance justification.

WHAT IT PRODUCES

Element copy inventory with genomic positions, structural integrity assessment per element, expression quantification, exogenous sequence discrimination, and burden quantification supporting clearance justification.

WHAT IT DETECTS

Endogenous retroviral element copies in the host genome with their integration positions, transcriptional expression of those elements, evidence of intact open reading frames, and the distinction between endogenous elements and exogenous infection.

WHAT IT DOES NOT DETECT

It does not establish particle infectivity, which requires biological assays. Whether an expressed element produces a competent particle cannot be determined from sequence and expression alone.

WHAT MAKES IT HARD

Retroviral elements are repetitive and highly similar, so mapping copies individually rather than as an aggregate requires careful handling of multi-mapping reads. Distinguishing an intact expressed element from a degenerate one determines whether it matters at all.

WHO DOES THIS JOB

Bioinformatician in viral safety, senior level, working with virology and regulatory affairs.

WHAT YOU CAN DO AFTERWARDS

Characterise endogenous retroviral elements and their expression, discriminate them from exogenous sequence, and support viral clearance justification.

WHO HIRES FOR IT

Biologics manufacturers using rodent cell lines, viral safety functions, contract testing organisations, and regulatory consultancies.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
18	88 hours · 15 days	Prior structural variant and expression analysis experience	2,05,000

COMPLETE WORKFLOW — 18 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Assessment scope confirmation	Element classes in scope and regulatory expectations confirmed	Protocol review	Assessment scope
02	Host genome reference preparation	Species-appropriate reference with element annotation prepared	Reference preparation	Prepared reference
03	Element reference assembly	Retroviral element reference sequences assembled for classification	Reference assembly	Element references
04	Data receipt and quality assessment	Sequence integrity and quality verified	md5sum, FastQC	QC acceptance record
05	Genome alignment	Reads aligned with multi-mapping handled explicitly given element similarity	bwa-mem2	Aligned BAM
06	Multi-mapping resolution	Reads mapping to multiple element copies resolved or assigned probabilistically	Multi-mapping methods	Resolved assignments
07	Element copy identification	Individual element copies identified with genomic positions	Copy identification	Element inventory
08	Element structure assessment	Each copy assessed for structural completeness and reading frame integrity	Structure analysis	Structural assessment
09	Intact element identification	Copies with intact reading frames identified as those of potential concern	Reading frame analysis	Intact element set
10	Expression data alignment	Transcriptomic data aligned with element-aware settings	STAR, element-aware alignment	Aligned RNA data
11	Element expression quantification	Expression quantified per element copy where mapping permits	Expression quantification	Expression profile
12	Aggregate expression assessment	Total element expression quantified where copies cannot be resolved	Aggregate quantification	Aggregate expression
13	Exogenous sequence discrimination	Exogenous retroviral sequence distinguished from endogenous elements	Discrimination analysis	Exogenous assessment
14	Particle-forming potential assessment	Structural and expression evidence combined to assess particle potential	Combined assessment	Potential assessment
15	Burden quantification	Element burden quantified for viral clearance justification	Burden calculation	Burden quantification
16	Infectivity limitation documentation	The inability of sequence data to establish infectivity documented	Limitation template	Infectivity caveat
17	Clearance justification support	Burden data assembled to support viral clearance validation	Support assembly	Clearance support data
18	Assessment report assembly	Inventory, structure, expression and burden assembled	Report assembly	Element assessment report

STEPS IN WORKFLOW

18

TRAINING DURATION

88 h · 15 d

PROGRAMME CODE

SBQ-07

TRAINING FEE (INR)

2,05,000

WHY IT IS DONE

Animal-derived and complex biological raw materials are the historical route by which adventitious agents entered manufacturing. Screening them before use is cheaper and safer than detecting the consequence in a bulk harvest.

WHAT IT ANSWERS

Does this raw material lot contain adventitious agent sequence, and does it qualify for use in manufacturing?

WHEN IT IS DONE

In qualification of new raw material suppliers, in lot-by-lot screening of high-risk materials, and in investigating a contamination traced to materials.

WHAT TRIGGERS IT

New supplier or material qualification, receipt of a high-risk material lot, or a contamination investigation implicating raw materials.

HOW IT IS DONE

Material matrix is characterised because it determines extraction and background, nucleic acid is extracted with controls processed identically, and reads are classified against comprehensive references. Source species is verified for animal-derived materials, detections are assessed against controls, and sensitivity is established by spike-in in the same matrix.

WHAT IT PRODUCES

Matrix characterisation, classification results with species assignment, source species verification, control comparison, matrix-specific spike-in sensitivity, and a material qualification recommendation.

WHAT IT DETECTS

Viral, bacterial, fungal and mycoplasmal sequence in the material, species identification where references permit, source species verification for animal-derived materials, and abundance relative to the material matrix.

WHAT IT DOES NOT DETECT

It detects nucleic acid rather than viable agent, and inactivated material may retain detectable sequence from a treated agent. Highly processed materials may carry degraded sequence that is undetectable while the original agent was present.

WHAT MAKES IT HARD

Sensitivity is matrix-dependent and cannot be transferred between material types, so spike-in must be performed in the actual matrix being tested. Source species verification frequently reveals undeclared substitution in animal-derived materials.

WHO DOES THIS JOB

Bioinformatician in raw material quality or biosafety, mid level.

WHAT YOU CAN DO AFTERWARDS

Screen raw materials for adventitious agents with matrix-appropriate sensitivity, verify source species, and support supplier qualification.

WHO HIRES FOR IT

Biologics manufacturers, raw material suppliers, contract testing organisations, and supply chain quality functions.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
17	82 hours · 14 days	Prior metagenomic analysis experience	1,90,000

COMPLETE WORKFLOW — 17 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Material scope and risk assessment	Material type, source and adventitious agent risk assessed	Risk assessment	<i>Material risk profile</i>
02	Testing protocol confirmation	Applicable testing scope and acceptance criteria confirmed	Protocol review	<i>Testing protocol</i>
03	Matrix characterisation	The material matrix characterised as it determines extraction and background	Matrix analysis	<i>Matrix characterisation</i>
04	Extraction method selection	Extraction selected and qualified for the specific material matrix	Method selection	<i>Selected extraction</i>
05	Matrix-specific spike preparation	Spike controls prepared in the actual material matrix	Spike preparation	<i>Matrix spike controls</i>
06	Sample and control extraction	Material and controls extracted in parallel under identical conditions	Extraction	<i>Prepared material</i>
07	Data receipt and quality assessment	Sequence integrity and quality verified per sample and control	md5sum, FastQC	<i>QC acceptance record</i>
08	Source species verification	Declared source species verified from the sequence data	Species analysis	<i>Species verification</i>
09	Species substitution assessment	Undeclared species substitution assessed for animal-derived materials	Substitution analysis	<i>Substitution assessment</i>
10	Host and matrix depletion	Expected source sequence depleted before agent classification	Depletion	<i>Depleted reads</i>
11	Agent classification	Remaining reads classified against comprehensive references	Classification	<i>Agent classification</i>
12	Reagent background subtraction	Detections assessed against reagent controls processed identically	Background analysis	<i>Corrected results</i>
13	Matrix-specific sensitivity determination	Sensitivity determined from spike recovery in the actual matrix	Recovery analysis	<i>Matrix sensitivity</i>
14	Detection threshold application	Detections assessed against the reporting threshold	Threshold criteria	<i>Threshold decisions</i>
15	Targeted confirmation	Detections confirmed by targeted analysis where above threshold	Confirmation analysis	<i>Confirmation results</i>
16	Qualification determination	A material qualification determination made against criteria	Determination criteria	<i>Qualification determination</i>
17	Screening report assembly	Species verification, detections and sensitivity assembled	Report assembly	<i>Screening report</i>

STEPS IN WORKFLOW

17

TRAINING DURATION

82 h · 14 d

PROGRAMME CODE

SBQ-08

TRAINING FEE (INR)

1,90,000

Transgene Integrity & Copy Number Analysis

(Transgene integration and copy number characterisation · Transgene characterisation)

WHY IT IS DONE

The transgene is the productive element of a recombinant cell line, and its copy number, integration structure and sequence integrity determine both productivity and product consistency. Characterising it is a defined expectation for cell substrate documentation.

WHAT IT ANSWERS

How many transgene copies does this cell line carry, where are they integrated, and is the transgene sequence intact and unarranged?

WHEN IT IS DONE

In master cell bank characterisation, in clone selection during cell line development, in stability studies, and in investigating productivity change.

WHAT TRIGGERS IT

Cell bank characterisation, clone selection, a stability study timepoint, or an unexplained productivity change.

HOW IT IS DONE

Copy number is quantified from coverage against the host background, integration junctions are recovered and mapped to host coordinates, and integration structure is reconstructed including tandem arrangements. Transgene sequence is compared against the intended construct across its full length, and flanking context is characterised.

WHAT IT PRODUCES

Copy number quantification, integration site coordinates with host gene context, integration structure reconstruction, full-length construct sequence verification, and a characterisation record for the substrate documentation.

WHAT IT DETECTS

Transgene copy number, integration site positions in the host genome, integration structure including concatemers and rearrangements, transgene sequence integrity against the construct, and flanking host sequence context.

WHAT IT DOES NOT DETECT

It does not measure expression level or productivity, which are separate attributes. Copy number alone does not predict productivity, since integration site context frequently matters more than copy count.

WHAT MAKES IT HARD

Tandem concatemeric integration is common and is easily reported as multiple independent sites unless junction structure is reconstructed. Copy number from coverage requires careful normalisation because the transgene is present at a very different abundance from the host genome.

WHO DOES THIS JOB

Bioinformatician in cell line development or quality, mid to senior level.

WHAT YOU CAN DO AFTERWARDS

Quantify transgene copy number, map and reconstruct integration structure, and verify construct sequence integrity.

WHO HIRES FOR IT

Biologics manufacturers, cell line development companies, contract development organisations, and quality functions.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
19	92 hours · 16 days	Prior structural variant and coverage analysis experience	2,10,000

COMPLETE WORKFLOW — 19 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Characterisation scope confirmation	Attributes to be characterised and their acceptance criteria confirmed	Protocol review	<i>Characterisation scope</i>
02	Construct reference registration	The intended construct sequence registered as the analytical reference	Sequence registry	<i>Registered construct</i>
03	Host reference preparation	Host genome reference prepared and version locked	Reference preparation	<i>Prepared reference</i>
04	Data receipt and quality assessment	Sequence integrity and quality verified against criteria	md5sum, FastQC	<i>QC acceptance record</i>
05	Combined reference alignment	Reads aligned against host and construct together	bwa-mem2, combined reference	<i>Aligned BAM</i>
06	Coverage normalisation	Construct coverage normalised against host genome background	Normalisation	<i>Normalised coverage</i>
07	Copy number quantification	Transgene copy number quantified with confidence bounds	Copy number quantification	<i>Copy number result</i>
08	Junction recovery	Vector to host junctions recovered from chimeric read evidence	Junction analysis	<i>Junction set</i>
09	Junction read-level review	Every candidate junction reviewed at read level to exclude artefact	IGV review	<i>Reviewed junctions</i>
10	Integration site mapping	Junctions mapped to host genomic coordinates	Coordinate mapping	<i>Integration sites</i>
11	Site count determination	The number of independent integration sites determined	Site analysis	<i>Site count</i>
12	Concatemer detection	Tandem concatemeric arrangements detected and distinguished from separate sites	Structure analysis	<i>Concatemer findings</i>
13	Integration structure reconstruction	Full integration architecture reconstructed per site	Structure reconstruction	<i>Integration structure</i>
14	Construct sequence verification	Integrated construct sequence verified across its full length	Sequence comparison	<i>Sequence verification</i>
15	Rearrangement and truncation assessment	Rearranged or truncated construct copies identified	Structural analysis	<i>Rearrangement findings</i>
16	Flanking context annotation	Host sequence and genes flanking each site annotated	Annotation	<i>Flanking context</i>
17	Site stability risk assessment	Integration site context assessed for stability risk	Risk assessment	<i>Site risk assessment</i>
18	Baseline archiving	The characterisation archived as the baseline for stability comparison	Archive procedure	<i>Archived baseline</i>
19	Characterisation report assembly	Copy number, sites, structure and sequence verification assembled	Report assembly	<i>Transgene characterisation report</i>

STEPS IN WORKFLOW

TRAINING DURATION

PROGRAMME CODE

TRAINING FEE (INR)

19

92 h · 16 d

SBQ-09

2,10,000

WHY IT IS DONE

Every recombinant product traces back to a construct sequence, and an undetected error in that sequence propagates into the product itself. Verification confirms that what was built matches what was designed before it becomes the basis of a manufacturing process.

WHAT IT ANSWERS

Does this expression construct match its intended design sequence exactly, and is every functional element present and correctly positioned?

WHEN IT IS DONE

At construct completion before cell line development, at plasmid bank establishment, and when a construct is transferred between sites or organisations.

WHAT TRIGGERS IT

Construct completion, plasmid bank establishment, a technology transfer, or a discrepancy investigation.

HOW IT IS DONE

The design sequence is registered as the reference, the construct is sequenced with a method resolving the full length including repetitive elements, and the assembled sequence is compared base by base against the design. Functional elements are verified for presence and boundaries, and any discrepancy is characterised and assessed for consequence.

WHAT IT PRODUCES

Full construct sequence assembly, base-level comparison against the design, functional element verification, discrepancy characterisation with consequence assessment, and a verification record for the construct history.

WHAT IT DETECTS

Full construct sequence with comparison against the design, sequence variants and their positions, functional element presence and boundaries, backbone integrity, and unintended sequence such as residual cloning scars.

WHAT IT DOES NOT DETECT

It does not establish that the construct functions as intended, which requires expression testing. Structural features such as supercoiling and methylation state are outside its scope.

WHAT MAKES IT HARD

Repetitive elements such as tandem promoters and terminators defeat short-read assembly and are exactly where errors occur, so the sequencing method must resolve them. Discrepancies require consequence assessment rather than a pass-fail decision, since many are silent.

WHO DOES THIS JOB

Bioinformatician in molecular biology or quality, mid level.

WHAT YOU CAN DO AFTERWARDS

Verify construct sequences against design at base resolution, resolve repetitive elements, and assess the consequence of discrepancies.

WHO HIRES FOR IT

Biologics and cell line development functions, gene therapy developers, plasmid manufacturers, and contract development organisations.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
17	82 hours · 14 days	Prior sequence analysis and assembly experience	1,90,000

COMPLETE WORKFLOW — 17 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Design sequence registration	The intended design sequence registered as the verification reference	Sequence registry	<i>Registered design</i>
02	Functional element annotation	Functional elements and their expected boundaries annotated on the design	Annotation	<i>Element map</i>
03	Sequencing method selection	A method resolving the full construct including repeats selected	Method selection	<i>Selected method</i>
04	Data receipt and quality assessment	Sequence integrity and run quality verified	md5sum, run QC	<i>QC acceptance record</i>
05	Read preprocessing	Reads preprocessed with settings appropriate to the method	Preprocessing	<i>Processed reads</i>
06	De novo assembly	The construct assembled independently of the design to avoid reference bias	Assembly	<i>Assembled construct</i>
07	Assembly quality assessment	Assembly completeness and consensus quality assessed	Assembly QC	<i>Assembly quality</i>
08	Circular structure resolution	Circular structure resolved with the origin defined consistently	Circularisation	<i>Resolved structure</i>
09	Design comparison	The assembly compared base by base against the design sequence	Sequence comparison	<i>Comparison results</i>
10	Variant identification	Sequence differences identified with positions and types recorded	Variant identification	<i>Variant list</i>
11	Repetitive region verification	Repetitive elements verified specifically since assembly errors concentrate there	Targeted verification	<i>Repeat verification</i>
12	Functional element verification	Each functional element verified for presence, sequence and boundaries	Element verification	<i>Element verification</i>
13	Backbone integrity verification	Backbone elements verified as intact and correctly positioned	Backbone analysis	<i>Backbone verification</i>
14	Unintended sequence detection	Cloning scars and unintended insertions detected	Sequence analysis	<i>Unintended sequence findings</i>
15	Consequence assessment	Each discrepancy assessed for functional consequence	Assessment framework	<i>Consequence assessment</i>
16	Verification determination	A verification determination made against acceptance criteria	Determination criteria	<i>Verification determination</i>
17	Verification record assembly	Assembly, comparison, discrepancies and determination assembled	Report assembly	<i>Verification record</i>

STEPS IN WORKFLOW

17

TRAINING DURATION

82 h · 14 d

PROGRAMME CODE

SBQ-10

TRAINING FEE (INR)

1,90,000

Producer Cell Line Clone Selection Analysis

(Clone selection genomic screening · Clone selection analysis)

WHY IT IS DONE

Selecting a production clone commits a product to that cell line for its commercial life, and the clone chosen on early productivity is frequently not the one that remains stable. Genomic screening at selection avoids discovering that years later.

WHAT IT ANSWERS

Which candidate clone has the integration structure, copy number and clonal purity most likely to remain stable through commercial manufacture?

WHEN IT IS DONE

During cell line development at the clone selection stage, and when re-selecting after a stability failure.

WHAT TRIGGERS IT

Completion of a cloning campaign requiring selection, or a stability failure requiring re-selection.

HOW IT IS DONE

Each candidate clone is characterised for copy number, integration structure and site context, and clonal purity is assessed for evidence that the population derives from a single cell. Integration site genomic environment is assessed for stability risk, clones are compared, and selection risk is ranked alongside productivity data.

WHAT IT PRODUCES

Per-clone copy number and integration characterisation, integration site context assessment, clonality evidence per clone, comparative stability risk ranking, and a selection recommendation with genomic rationale.

WHAT IT DETECTS

Transgene copy number and integration structure per clone, integration site genomic context including transcriptional environment, clonal purity and evidence of monoclonality, and early indicators of instability.

WHAT IT DOES NOT DETECT

It does not measure productivity or product quality, which drive selection alongside genomic criteria. Long-term stability cannot be predicted with certainty from early genomic characterisation, only risk-ranked.

WHAT MAKES IT HARD

Monoclonality evidence is a regulatory expectation and genomic data contribute to it, but no single method proves it, so the argument is built from converging evidence. Integration site context predicts stability better than copy number and is routinely ignored in favour of the simpler metric.

WHO DOES THIS JOB

Bioinformatician in cell line development, mid to senior level, working with process development.

WHAT YOU CAN DO AFTERWARDS

Characterise clones genomically at selection, assess integration site risk, contribute to monoclonality arguments, and rank clones on stability risk.

WHO HIRES FOR IT

Biologics manufacturers, cell line development companies, contract development organisations, and cell engineering functions.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
20	96 hours · 16 days	Prior copy number and clonality analysis experience	2,20,000

COMPLETE WORKFLOW — 20 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Selection criteria definition	Genomic selection criteria defined alongside productivity criteria	Criteria specification	Selection criteria
02	Clone inventory and tracking	Candidate clones inventoried with cloning history and tracking verified	Inventory, tracking	Clone inventory
03	Data receipt and quality assessment	Sequence integrity and quality verified per clone	md5sum, FastQC	QC acceptance record
04	Alignment and processing	All clones aligned and processed identically	bwa-mem2, Picard	Analysis-ready BAMs
05	Coverage consistency assessment	Coverage confirmed comparable across clones for valid comparison	mosdepth	Coverage consistency
06	Copy number quantification	Transgene copy number quantified per clone	Copy number quantification	Per-clone copy number
07	Junction recovery per clone	Integration junctions recovered for every candidate clone	Junction analysis	Per-clone junctions
08	Integration site mapping	Sites mapped to host coordinates for each clone	Coordinate mapping	Per-clone sites
09	Integration structure reconstruction	Integration architecture reconstructed per clone	Structure reconstruction	Per-clone structures
10	Site genomic context annotation	Genomic and transcriptional context of each site annotated	Annotation	Context annotation
11	Transcriptional environment assessment	Site chromatin and transcriptional environment assessed for expression stability	Environment analysis	Environment assessment
12	Fragile and repetitive site flagging	Sites in repetitive or unstable genomic regions flagged	Region analysis	Site risk flags
13	Clonal purity assessment	Each clone assessed for evidence of a single-cell origin	Clonality analysis	Clonal purity
14	Heterogeneity detection	Subpopulation heterogeneity within each clone detected	Heterogeneity analysis	Heterogeneity findings
15	Monoclonality evidence assembly	Genomic evidence assembled to contribute to monoclonality assurance	Evidence assembly	Monoclonality evidence
16	Construct integrity verification	Construct sequence integrity verified in each clone	Sequence comparison	Integrity verification
17	Early instability indicator assessment	Early indicators of instability assessed across clones	Indicator analysis	Instability indicators
18	Comparative risk ranking	Clones ranked on combined genomic stability risk	Ranking framework	Risk ranking
19	Selection recommendation	A selection recommendation made integrating genomic and productivity evidence	Recommendation framework	Selection recommendation
20	Selection report assembly	Per-clone characterisation, evidence and ranking assembled	Report assembly	Clone selection report

STEPS IN WORKFLOW	TRAINING DURATION	PROGRAMME CODE	TRAINING FEE (INR)
20	96 h · 16 d	SBQ-11	2,20,000

Host Cell DNA Residual Analysis

(Residual host cell DNA characterisation · Residual DNA analysis)

WHY IT IS DONE

Residual host cell DNA in a purified biologic is a defined impurity with limits on both quantity and fragment size, because the risk it poses depends on whether intact functional sequences survive. Sequencing characterises what quantitative assays only measure in bulk.

WHAT IT ANSWERS

What host cell DNA remains in this purified product, what is its size distribution, and does any functionally relevant sequence survive intact?

WHEN IT IS DONE

In process characterisation to demonstrate clearance, in supporting residual DNA specifications, and in investigating a residual DNA result that exceeds expectation.

WHAT TRIGGERS IT

Process characterisation, a specification justification requirement, or an out-of-expectation residual DNA measurement.

HOW IT IS DONE

Residual DNA is recovered from purified product with a method preserving fragment size, sequenced with size information retained, and fragments are mapped to the host genome. Size distribution is characterised, genomic origin is analysed for regional enrichment, and functionally relevant sequences are assessed for intactness.

WHAT IT PRODUCES

Fragment size distribution, genomic origin mapping with regional enrichment, intactness assessment for functionally relevant sequences, construct sequence detection, and a characterisation supporting risk assessment.

WHAT IT DETECTS

Host DNA fragment size distribution, genomic origin of surviving fragments, whether oncogenic or coding sequences survive intact, enrichment of particular genomic regions, and vector or construct sequence carried alongside.

WHAT IT DOES NOT DETECT

It characterises what is present rather than quantifying total mass, which requires an orthogonal quantitative method. It also cannot establish biological risk, only the sequence characteristics that risk assessment depends on.

WHAT MAKES IT HARD

Recovery from purified product is technically demanding because the DNA is present at very low mass among abundant protein, and any recovery bias distorts the size distribution that the whole analysis depends on. Fragment size is the risk-relevant attribute and is easily lost through library preparation.

WHO DOES THIS JOB

Bioinformatician in analytical development or quality control, mid to senior level.

WHAT YOU CAN DO AFTERWARDS

Characterise residual host DNA including fragment size and genomic origin, assess intactness of relevant sequences, and support risk-based specifications.

WHO HIRES FOR IT

Biologics manufacturers, analytical development functions, contract testing laboratories, and regulatory affairs groups.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
18	86 hours · 15 days	Prior sequencing and quantification analysis experience	2,00,000

Host Cell DNA Residual Analysis

(Residual host cell DNA characterisation · Residual DNA analysis)

COMPLETE WORKFLOW — 18 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Analysis scope confirmation	Attributes required and their regulatory basis confirmed	Protocol review	Analysis scope
02	Recovery method qualification	The recovery method qualified for the product matrix with bias assessed	Method qualification	Qualified method
03	Size preservation verification	Library preparation confirmed to preserve fragment size information	Method verification	Size preservation record
04	Sample and control preparation	Product and process controls prepared under identical conditions	Preparation	Prepared material
05	Data receipt and quality assessment	Sequence integrity and quality verified	md5sum, FastQC	QC acceptance record
06	Host genome alignment	Recovered fragments aligned to the host cell genome	bwa-mem2	Aligned reads
07	Non-host classification	Reads not of host origin classified and characterised separately	Classification	Non-host classification
08	Fragment size distribution analysis	Fragment length distribution characterised across the recovered population	Length analysis	Size distribution
09	Size class quantification	Fragments quantified within regulatory-relevant size classes	Class quantification	Size class proportions
10	Genomic origin mapping	Fragment origin mapped across the host genome	Coordinate mapping	Origin map
11	Regional enrichment analysis	Enrichment of particular genomic regions assessed statistically	Enrichment analysis	Enrichment findings
12	Coding sequence intactness assessment	Whether coding sequences survive intact assessed	Intactness analysis	Coding intactness
13	Oncogenic sequence assessment	Survival of sequences of specific regulatory concern assessed	Targeted analysis	Oncogenic sequence findings
14	Construct sequence detection	Vector or construct sequence carried with host DNA detected	Sequence detection	Construct findings
15	Recovery bias assessment	Method recovery bias assessed against spiked size standards	Bias analysis	Bias assessment
16	Quantitative method reconciliation	Sequencing results reconciled against orthogonal quantitative measurement	Reconciliation	Reconciliation record
17	Risk-relevant characterisation summary	Attributes relevant to risk assessment summarised	Summary assembly	Risk characterisation
18	Characterisation report assembly	Size distribution, origin, intactness and reconciliation assembled	Report assembly	Residual DNA report

STEPS IN WORKFLOW	TRAINING DURATION	PROGRAMME CODE	TRAINING FEE (INR)
18	86 h · 15 d	SBQ-12	2,00,000

Viral Vector Genome Integrity Analysis

(Vector genome integrity and impurity characterisation · Vector integrity QC)

WHY IT IS DONE

A viral vector preparation is a mixture of intended genomes and a range of incomplete, rearranged and foreign species. The proportion of full-length intended genome is a direct determinant of potency, and the impurities are a direct safety consideration.

WHAT IT ANSWERS

What proportion of packaged genomes in this vector lot are complete and correct, and what impurity species make up the remainder?

WHEN IT IS DONE

In vector batch release, in process development and comparability, in investigating potency variability, and in supporting regulatory submissions.

WHAT TRIGGERS IT

A vector batch requiring release, a process change, a potency deviation, or a submission requirement.

HOW IT IS DONE

Unencapsidated DNA is removed before extraction so that packaged content is measured, sequencing is performed with a method resolving full-length genomes and terminal repeats, and reads are classified against vector, host and helper references. Genome integrity is assessed across the full construct, truncation positions are mapped, and impurity proportions are quantified against specification.

WHAT IT PRODUCES

Packaged content classification, full-length genome proportion, truncation position mapping, terminal repeat integrity assessment, reverse packaged and recombinant species quantification, and specification compliance.

WHAT IT DETECTS

Full-length genome proportion, truncated species with termination positions mapped, terminal repeat integrity, reverse packaged and inverted species, recombinant genomes, and packaged host and helper sequence.

WHAT IT DOES NOT DETECT

It does not measure infectious titre, capsid content ratio or transduction efficiency, which require separate methods. Encapsidated and surface-associated DNA are indistinguishable without appropriate sample treatment before extraction.

WHAT MAKES IT HARD

Terminal repeats are structurally difficult to sequence and are precisely where integrity failures occur, so a method that cannot resolve them is not characterising the attribute that matters. Skipping nuclease treatment produces impurity figures that describe the buffer rather than the vector.

WHO DOES THIS JOB

Bioinformatician in vector analytical development, senior level, working with process development and quality control.

WHAT YOU CAN DO AFTERWARDS

Characterise packaged vector genomes, resolve terminal repeats, quantify impurity species, and support release and comparability decisions.

WHO HIRES FOR IT

Gene therapy developers, vector contract manufacturers, analytical development functions, and quality control laboratories.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
22	106 hours · 18 days	Prior sequencing and assembly experience	2,45,000

COMPLETE WORKFLOW — 22 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Construct and specification registration	Vector genome sequence and release specifications registered	Registry	<i>Registered references</i>
02	Reference set assembly	Vector, host and helper references assembled for classification	Reference assembly	<i>Reference set</i>
03	Nuclease treatment verification	Removal of unencapsidated DNA verified before extraction	Treatment verification	<i>Treatment record</i>
04	Extraction and input quantification	Packaged nucleic acid extracted and quantified	Extraction, quantification	<i>Input record</i>
05	Sequencing method selection	A method resolving full-length genomes and terminal repeats selected	Method selection	<i>Selected method</i>
06	Data receipt and quality assessment	Sequence integrity and run quality verified	md5sum, run QC	<i>QC acceptance record</i>
07	Read classification	Reads classified against vector, host and helper references	Classification	<i>Classified reads</i>
08	Vector genome alignment	Vector-classified reads aligned to the construct reference	Alignment	<i>Aligned reads</i>
09	Coverage profile analysis	Coverage analysed across the construct for uniformity and gaps	Coverage analysis	<i>Coverage profile</i>
10	Full-length species quantification	Complete genome species quantified as a proportion of packaged content	Length analysis	<i>Full-length proportion</i>
11	Truncation position mapping	Termination positions of truncated species mapped	Truncation analysis	<i>Truncation map</i>
12	Truncated species quantification	Truncated species quantified by size class	Quantification	<i>Truncation profile</i>
13	Terminal repeat analysis	Terminal repeat structure and integrity analysed specifically	Specialised analysis	<i>Terminal repeat assessment</i>
14	Terminal repeat aberration detection	Deletions and rearrangements within terminal repeats detected	Aberration analysis	<i>Aberration findings</i>
15	Reverse packaged species assessment	Reverse packaged sequence identified and quantified	Orientation analysis	<i>Reverse packaging findings</i>
16	Recombinant genome detection	Recombinant and rearranged genome species detected	Recombination analysis	<i>Recombinant findings</i>
17	Host DNA packaging quantification	Packaged host sequence quantified with genomic origin characterised	Host analysis	<i>Host packaging profile</i>
18	Helper sequence quantification	Helper plasmid packaging quantified	Helper analysis	<i>Helper profile</i>
19	Sequence variant analysis	Variants within the vector genome called against error background	Variant analysis	<i>Variant profile</i>
20	Specification comparison	All measured attributes compared against release specifications	Specification comparison	<i>Compliance assessment</i>
21	Batch trend comparison	Results compared against historical batches for process consistency	Trend analysis	<i>Trend comparison</i>
22	Release report assembly	Integrity, impurities, compliance and trends assembled	Report assembly	<i>Vector integrity report</i>

STEPS IN WORKFLOW

22

TRAINING DURATION

106 h · 18 d

PROGRAMME CODE

SBQ-13

TRAINING FEE (INR)

2,45,000

Capsid & Packaging Composition Analysis

(Vector packaging composition characterisation · Packaging composition analysis)

WHY IT IS DONE

What a vector particle contains determines both how much product reaches the target and what unintended material is delivered with it. Characterising packaging composition is how a manufacturer explains potency variation between lots that look identical by titre.

WHAT IT ANSWERS

What proportion of packaged material is intended vector genome, and what proportion is host, helper or other sequence?

WHEN IT IS DONE

In process development to reduce impurity packaging, in batch comparability, in potency investigation, and in supporting specification setting.

WHAT TRIGGERS IT

A process change affecting packaging, a potency investigation, a comparability exercise, or specification development.

HOW IT IS DONE

Unencapsidated material is removed, packaged nucleic acid is extracted with size preserved, and reads are classified against all candidate source references. Proportions are quantified with normalisation for genome size, host packaging is analysed for regional origin, and composition is compared across lots and process conditions.

WHAT IT PRODUCES

Packaged content proportions by source with size normalisation, host DNA regional origin distribution, packaged species size distribution, process condition comparison, and specification support data.

WHAT IT DETECTS

Relative proportions of vector, host cell, helper plasmid and other packaged sequence, genomic origin distribution of host DNA packaged, size distribution of packaged species, and process-related differences in composition.

WHAT IT DOES NOT DETECT

It does not measure empty to full particle ratio, which requires biophysical methods, and it does not establish which packaged species are biologically consequential.

WHAT MAKES IT HARD

Proportions must be normalised for the very different genome sizes of vector, helper and host or the composition figures are meaningless. Recovery bias during extraction distorts the size distribution that composition analysis depends on.

WHO DOES THIS JOB

Bioinformatician in vector analytical development, mid to senior level.

WHAT YOU CAN DO AFTERWARDS

Quantify packaging composition with correct normalisation, characterise host DNA packaging, and support process optimisation and specification setting.

WHO HIRES FOR IT

Gene therapy developers, vector manufacturers, analytical development functions, and process development groups.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
20	98 hours · 17 days	Prior vector analysis experience	2,25,000

Capsid & Packaging Composition Analysis

(Vector packaging composition characterisation · Packaging composition analysis)

COMPLETE WORKFLOW — 20 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Analysis scope confirmation	Packaging attributes in scope and their intended use confirmed	Protocol review	<i>Analysis scope</i>
02	Reference set assembly	Vector, host, helper and other candidate references assembled	Reference assembly	<i>Reference set</i>
03	Nuclease treatment verification	Unencapsidated material removal verified before extraction	Treatment verification	<i>Treatment record</i>
04	Size-preserving extraction	Nucleic acid extracted with fragment size information preserved	Extraction	<i>Extracted material</i>
05	Data receipt and quality assessment	Sequence integrity and run quality verified	md5sum, run QC	<i>QC acceptance record</i>
06	Read classification	Reads assigned to vector, host, helper or unassigned categories	Classification	<i>Classified reads</i>
07	Ambiguous read handling	Reads mapping to more than one source handled under a documented rule	Ambiguity handling	<i>Ambiguity record</i>
08	Genome size normalisation	Proportions normalised for the differing genome sizes of each source	Normalisation	<i>Normalised proportions</i>
09	Packaging proportion quantification	Proportions of each packaged species quantified with bounds	Quantification	<i>Packaging proportions</i>
10	Host DNA origin mapping	Packaged host sequence mapped across the host genome	Coordinate mapping	<i>Host origin map</i>
11	Host regional enrichment analysis	Enrichment of particular host regions in packaged DNA assessed	Enrichment analysis	<i>Regional enrichment</i>
12	Packaged size distribution analysis	Size distribution of packaged species analysed by source	Length analysis	<i>Size distribution</i>
13	Helper sequence characterisation	Packaged helper sequence characterised by region and size	Helper analysis	<i>Helper characterisation</i>
14	Unassigned sequence characterisation	Sequence matching no reference assembled and characterised	Assembly, characterisation	<i>Unassigned findings</i>
15	Recovery bias assessment	Extraction and library bias assessed against size standards	Bias analysis	<i>Bias assessment</i>
16	Process condition comparison	Composition compared across process conditions or batches	Comparative analysis	<i>Condition comparison</i>
17	Historical trend comparison	Results compared against historical composition data	Trend analysis	<i>Trend comparison</i>
18	Specification support analysis	Data assembled to support composition specification setting	Support analysis	<i>Specification support</i>
19	Method limitation documentation	Attributes outside the method's scope documented explicitly	Limitation template	<i>Limitation statement</i>
20	Composition report assembly	Proportions, origin, distribution and comparisons assembled	Report assembly	<i>Packaging composition report</i>

STEPS IN WORKFLOW

20

TRAINING DURATION

98 h · 17 d

PROGRAMME CODE

SBQ-14

TRAINING FEE (INR)

2,25,000

Plasmid Identity & Quality Analysis

(Plasmid sequence identity and quality testing · Plasmid QC)

WHY IT IS DONE

Plasmids are starting materials for vectors, mRNA and cell therapies, and an error or impurity in the plasmid propagates into everything made from it. Sequence-level quality testing is a starting material control with downstream consequences at every stage.

WHAT IT ANSWERS

Does this plasmid lot match its reference sequence, and what sequence-level impurities are present?

WHEN IT IS DONE

In plasmid lot release, in plasmid bank establishment, in supplier qualification, and in investigating a downstream deviation traced to plasmid material.

WHAT TRIGGERS IT

Plasmid lot receipt or release, bank establishment, supplier qualification, or a downstream deviation investigation.

HOW IT IS DONE

The reference sequence is registered, the plasmid is sequenced with a method resolving the full circular sequence including repetitive regions, and the assembly is compared base by base against the reference. Host genomic contamination is quantified, other plasmid species are detected, and variants are assessed for consequence.

WHAT IT PRODUCES

Full sequence assembly with reference comparison, variant identification with consequence assessment, structural integrity confirmation, host DNA contamination quantification, foreign plasmid detection, and release recommendation.

WHAT IT DETECTS

Full plasmid sequence with comparison against the reference, sequence variants and their positions, structural rearrangements and deletions, host genomic DNA contamination, and other plasmid species present.

WHAT IT DOES NOT DETECT

It does not measure supercoiled fraction, endotoxin or residual host protein, which are separate release attributes. Plasmid topology is outside sequencing scope entirely.

WHAT MAKES IT HARD

Circular assembly across repetitive backbone elements requires appropriate read lengths, and short-read assembly frequently produces a plausible but incorrect structure. Low-level variant populations within a plasmid preparation are real and require the sequencing error background to be characterised.

WHO DOES THIS JOB

Bioinformatician in quality control or analytical development, mid level.

WHAT YOU CAN DO AFTERWARDS

Verify plasmid sequence identity at base resolution, resolve circular and repetitive structure, and quantify sequence-level impurities.

WHO HIRES FOR IT

Plasmid manufacturers, gene therapy and mRNA developers, contract manufacturers, and quality control laboratories.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
17	82 hours · 14 days	Prior sequence assembly experience	1,90,000

COMPLETE WORKFLOW — 17 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Reference sequence registration	The reference plasmid sequence registered for comparison	Sequence registry	Registered reference
02	Specification confirmation	Applicable release specifications and attributes confirmed	Specification review	Specification set
03	Sequencing method selection	A method resolving the full circular sequence including repeats selected	Method selection	Selected method
04	Data receipt and quality assessment	Sequence integrity and run quality verified	md5sum, run QC	QC acceptance record
05	Read preprocessing	Reads preprocessed with settings appropriate to the method	Preprocessing	Processed reads
06	De novo assembly	The plasmid assembled independently of the reference	Assembly	Assembled plasmid
07	Circular structure resolution	Circularity resolved with the origin defined consistently	Circularisation	Resolved structure
08	Assembly quality assessment	Assembly completeness and consensus accuracy assessed	Assembly QC	Assembly quality
09	Reference comparison	Assembly compared base by base against the registered reference	Sequence comparison	Comparison results
10	Variant identification and consequence assessment	Differences identified and assessed for functional consequence	Variant analysis	Variant assessment
11	Structural integrity verification	Deletions, insertions and rearrangements assessed against the reference	Structural analysis	Structural verification
12	Repetitive region verification	Repetitive backbone elements verified specifically	Targeted verification	Repeat verification
13	Low-frequency variant analysis	Subpopulation variants assessed against the sequencing error background	Variant analysis	Subpopulation findings
14	Host genomic DNA quantification	Host cell genomic DNA contamination quantified	Contamination quantification	Host DNA level
15	Foreign plasmid detection	Other plasmid species present in the preparation detected	Sequence classification	Foreign plasmid findings
16	Specification comparison	All measured attributes compared against release specifications	Specification comparison	Compliance assessment
17	Release report assembly	Sequence verification, impurities and compliance assembled	Report assembly	Plasmid QC report

STEPS IN WORKFLOW

17

TRAINING DURATION

82 h · 14 d

PROGRAMME CODE

SBQ-15

TRAINING FEE (INR)

1,90,000

WHY IT IS DONE

For an mRNA product the sequence is the active substance, so identity confirmation and impurity characterisation are direct quality attributes rather than supporting information. Truncated and fragmented species affect both potency and immunogenicity.

WHAT IT ANSWERS

Does this mRNA lot match its intended sequence, and what impurity species are present at what proportions?

WHEN IT IS DONE

In batch release, in process development and comparability, in stability studies, and in investigating a potency or immunogenicity signal.

WHAT TRIGGERS IT

A batch requiring release, a process change, a stability timepoint, or a potency investigation.

HOW IT IS DONE

Sequencing chemistry is selected to resolve full-length species rather than fragment them, reads are aligned to the intended sequence, and identity is confirmed across the full length. Length distributions quantify truncation and fragmentation, variants are called against a characterised error background, tail length is profiled, and residual template is quantified.

WHAT IT PRODUCES

Full-length identity confirmation, truncation and fragmentation profiles with positions, variant analysis against error background, tail length distribution, residual template quantification, and specification comparison.

WHAT IT DETECTS

Full-length sequence identity, truncated species with termination positions, fragmented species proportions, sequence variants against the error background, poly-adenosine tail length distribution, and residual template sequence.

WHAT IT DOES NOT DETECT

It does not measure capping efficiency, modified nucleoside incorporation, secondary structure or lipid formulation attributes, all of which require separate methods. Sequencing chemistry introduces artefacts that mimic product impurities.

WHAT MAKES IT HARD

Separating genuine product impurity from sequencing artefact requires a characterised control and an explicit error model, because both produce identical signatures. Fragmenting library preparation destroys the length information the impurity analysis depends on.

WHO DOES THIS JOB

Bioinformatician in analytical development or quality control, mid to senior level.

WHAT YOU CAN DO AFTERWARDS

Confirm mRNA identity, quantify impurity species against error background, characterise tail length, and support release decisions.

WHO HIRES FOR IT

mRNA therapeutic and vaccine developers, contract manufacturers, analytical development functions, and quality control laboratories.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
20	96 hours · 16 days	Prior sequencing analysis and long-read basics	2,20,000

COMPLETE WORKFLOW — 20 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Construct sequence registration	The intended product sequence registered as the analytical reference	Sequence registry	Registered construct
02	Quality attribute confirmation	Sequence-related quality attributes and specifications confirmed	Specification review	Attribute set
03	Sequencing method selection	A method resolving full-length species selected	Method selection	Selected method
04	Control material preparation	Reference material processed identically to characterise method background	Control preparation	Control data
05	Data receipt and quality assessment	Sequence integrity and run quality verified	md5sum, run QC	QC acceptance record
06	Read preprocessing	Reads preprocessed with length information preserved	Preprocessing	Processed reads
07	Alignment to intended sequence	Reads aligned to the registered construct sequence	Alignment	Aligned reads
08	Coverage profile analysis	Coverage analysed across the full construct for gaps and bias	Coverage analysis	Coverage profile
09	Full-length identity confirmation	Identity confirmed across the entire construct length	Identity analysis	Identity confirmation
10	Read length distribution analysis	Length distribution analysed to characterise product species	Length analysis	Length distribution
11	Truncation position mapping	Termination positions of truncated species mapped	Truncation analysis	Truncation map
12	Truncated species quantification	Truncated species quantified as a proportion of total	Quantification	Truncation profile
13	Fragmentation quantification	Fragmented species quantified with size classes defined	Fragmentation analysis	Fragmentation profile
14	Method error background characterisation	Sequencing error background characterised from control material	Background analysis	Error background
15	Sequence variant analysis	Variants called against the characterised error background	Variant analysis	Variant profile
16	Poly-adenosine tail characterisation	Tail length distribution characterised where applicable	Tail analysis	Tail profile
17	Residual template quantification	Residual DNA template or plasmid sequence quantified	Template analysis	Residual template
18	Impurity classification	All impurity species classified by type and quantified	Classification	Impurity profile
19	Specification comparison	Measured attributes compared against specifications	Specification comparison	Compliance assessment
20	Release report assembly	Identity, impurities and compliance assembled	Report assembly	mRNA QC report

STEPS IN WORKFLOW	TRAINING DURATION	PROGRAMME CODE	TRAINING FEE (INR)
20	96 h · 16 d	SBQ-16	2,20,000

Cell Therapy Product Release Analysis

(Cell therapy genomic release testing · Cell therapy release testing)

WHY IT IS DONE

A cell therapy product is manufactured per patient and released against specification within days, which compresses characterisation into a window no other biologic faces. Genomic release attributes must be measurable, reproducible and fast.

WHAT IT ANSWERS

Does this cell therapy product meet its genomic release specifications, and is it safe and consistent enough to infuse?

WHEN IT IS DONE

At batch release for every patient product, in process comparability, and in investigating an out-of-specification release result.

WHAT TRIGGERS IT

A patient product requiring release, a process comparability exercise, or an out-of-specification investigation.

HOW IT IS DONE

Release attributes and their specifications are confirmed, patient identity concordance is verified between starting material and product, and vector copy number and transduction are quantified against specification. Adventitious agent and replication competent vector screening are performed within the release window, and composition attributes are measured with results assembled against specification.

WHAT IT PRODUCES

Patient identity concordance, vector copy number and transduced proportion against specification, adventitious agent screening results, replication competent vector assessment, composition attributes, and a release recommendation.

WHAT IT DETECTS

Vector copy number per cell, transduced cell proportion, product identity against the patient, adventitious agent status, replication competent vector signal where applicable, and product composition attributes.

WHAT IT DOES NOT DETECT

It does not measure potency or cytotoxic function, which require functional assays. The release window limits characterisation depth, so some attributes are assessed retrospectively rather than before infusion.

WHAT MAKES IT HARD

Turnaround is the binding constraint, since a result arriving after the infusion window has no release value. Patient identity concordance between apheresis material and final product is the control against the most serious possible error and must be built into the workflow rather than added.

WHO DOES THIS JOB

Bioinformatician in cell therapy quality control, senior level, operating within a release timeline.

WHAT YOU CAN DO AFTERWARDS

Operate genomic release testing for cell therapy within a compressed window, verify patient identity concordance, and make specification-based release recommendations.

WHO HIRES FOR IT

Cell therapy manufacturers, hospital-based manufacturing units, contract manufacturing organisations, and quality control functions.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
22	108 hours · 18 days	Prior single-cell and vector analysis experience	2,45,000

COMPLETE WORKFLOW — 22 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Release specification confirmation	Genomic release attributes and their specifications confirmed	Specification review	Release specifications
02	Release window recording	The available time from sampling to infusion recorded as the operating constraint	Clinical protocol	Release window
03	Starting material identity capture	Patient identity captured from apheresis starting material	Identity analysis	Starting material identity
04	Product sample custody verification	Custody from product to analytical sample verified	Custody procedure	Custody record
05	Data receipt and quality assessment	Sequence integrity and quality verified within the release window	md5sum, FastQC	QC acceptance record
06	Patient identity concordance	Product identity confirmed to match the starting material and patient	Concordance analysis	Identity concordance
07	Vector sequence verification	The transgene sequence verified against the intended construct	Sequence verification	Construct verification
08	Vector copy number quantification	Vector copy number per cell quantified against specification	Copy number quantification	Copy number result
09	Transduced proportion determination	Proportion of transduced cells determined	Transduction analysis	Transduced proportion
10	Integration site assessment	Integration site distribution assessed where the specification requires it	Integration analysis	Integration assessment
11	Clonal dominance assessment	Dominant integration clones assessed against safety criteria	Dominance analysis	Dominance assessment
12	Replication competent vector screening	Replication competent vector sequence screened for where applicable	Screening analysis	RCV screening result
13	Adventitious agent screening	Rapid adventitious agent screening performed within the window	Screening analysis	Agent screening result
14	Microbial detection	Bacterial and mycoplasmal sequence screened for within the window	Detection analysis	Microbial result
15	Control comparison	All detections assessed against controls processed in parallel	Control analysis	Control comparison
16	Product composition assessment	Cell composition attributes measured where specified	Composition analysis	Composition result
17	Specification comparison	Every attribute compared against its release specification	Specification comparison	Compliance assessment
18	Out-of-specification flagging	Any attribute outside specification flagged for immediate investigation	Flagging procedure	OOS flags
19	Turnaround verification	Result delivery verified against the release window commitment	Timing audit	Turnaround record
20	Release recommendation	A release recommendation made against the specification set	Recommendation framework	Release recommendation
21	Retrospective characterisation planning	Attributes for retrospective analysis after infusion identified	Planning procedure	Retrospective plan
22	Release report assembly	Identity, attributes, compliance and recommendation assembled	Report assembly	Release report

STEPS IN WORKFLOW

22

TRAINING DURATION

108 h · 18 d

PROGRAMME CODE

SBQ-17

TRAINING FEE (INR)

2,45,000

Microbial Identification for Contamination Investigation

(Contamination isolate identification and source tracing · Contamination investigation)

WHY IT IS DONE

When a contamination event occurs, identifying the organism is only the first question — the decisive one is where it came from. Whole genome comparison against environmental and historical isolates traces the route in a way that species identification alone never can.

WHAT IT ANSWERS

What organism caused this contamination, and does its genome match any isolate previously recovered from the facility, personnel or materials?

WHEN IT IS DONE

Immediately following a contamination event, in root cause investigation, and in demonstrating that corrective actions eliminated the source.

WHAT TRIGGERS IT

A sterility failure, a bioburden excursion, an environmental monitoring alert, or a media fill failure.

HOW IT IS DONE

The isolate is sequenced and assembled, species and strain are identified, and the genome is compared against the facility's environmental monitoring isolate archive. Relatedness is quantified at core genome level with thresholds interpreted against organism-specific diversity, persistence across time is assessed, and resistance determinants are characterised.

WHAT IT PRODUCES

Species and strain identification, genomic comparison against the isolate archive with relatedness quantified, persistence assessment, resistance determinant characterisation, source hypothesis with evidence, and corrective action support.

WHAT IT DETECTS

Species and strain identification of the isolate, genomic relatedness to environmental monitoring isolates, relatedness to historical facility isolates, antimicrobial and biocide resistance determinants, and persistence of a strain across time.

WHAT IT DOES NOT DETECT

It cannot establish the transmission direction between two related isolates, only that they are related. Absence of a match means the source was not sampled rather than that no source exists.

WHAT MAKES IT HARD

The comparison is only as good as the archive, so facilities without a systematic environmental isolate collection cannot trace anything. Relatedness thresholds differ by organism, and applying a single cutoff produces both false links and missed ones.

WHO DOES THIS JOB

Bioinformatician in microbiology or quality investigation, mid to senior level, working with microbiology and quality assurance.

WHAT YOU CAN DO AFTERWARDS

Identify contamination isolates to strain level, compare against facility archives, quantify relatedness correctly, and construct evidence-based source hypotheses.

WHO HIRES FOR IT

Biologics and sterile manufacturers, contract manufacturing organisations, quality microbiology functions, and pharmaceutical microbiology consultancies.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
18	86 hours · 15 days	Prior microbial genomics experience	2,00,000

Microbial Identification for Contamination Investigation

(Contamination isolate identification and source tracing · Contamination investigation)

COMPLETE WORKFLOW — 18 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Event scope and timeline capture	The contamination event, affected material and timeline captured	Event record	<i>Event record</i>
02	Isolate custody verification	Custody of the isolate from recovery to sequencing verified	Custody procedure	<i>Custody record</i>
03	Archive retrieval	Environmental and historical facility isolate genomes retrieved for comparison	Archive retrieval	<i>Comparison archive</i>
04	Data receipt and quality assessment	Sequence integrity and quality verified for the isolate	md5sum, FastQC	<i>QC acceptance record</i>
05	Genome assembly	The isolate genome assembled and quality assessed	Assembly, assembly QC	<i>Assembled genome</i>
06	Species identification	Species identified from whole genome comparison	Taxonomic analysis	<i>Species identification</i>
07	Strain typing	Strain assigned by the appropriate typing scheme for the organism	Typing analysis	<i>Strain assignment</i>
08	Core genome comparison	Core genome comparison performed against archive isolates	Core genome analysis	<i>Core comparison</i>
09	Relatedness quantification	Genomic distance quantified against archive isolates	Distance analysis	<i>Relatedness metrics</i>
10	Organism-specific threshold application	Relatedness assessed against thresholds appropriate to the organism's diversity	Threshold criteria	<i>Relatedness determination</i>
11	Phylogenetic placement	The isolate placed phylogenetically among archive and reference isolates	Phylogenetic analysis	<i>Phylogenetic placement</i>
12	Persistence assessment	Evidence of the strain persisting across time in the facility assessed	Temporal analysis	<i>Persistence assessment</i>
13	Resistance determinant characterisation	Antimicrobial and biocide resistance determinants characterised	Resistance analysis	<i>Resistance profile</i>
14	Biofilm and survival trait assessment	Traits supporting persistence in the facility environment assessed	Trait analysis	<i>Survival trait findings</i>
15	Source hypothesis construction	A source hypothesis constructed from the genomic and facility evidence	Hypothesis framework	<i>Source hypothesis</i>
16	Archive gap documentation	Locations and materials not represented in the archive documented	Gap analysis	<i>Archive gap record</i>
17	Corrective action support	Genomic evidence assembled to support corrective and preventive action	Support assembly	<i>Corrective action support</i>
18	Investigation report assembly	Identification, relatedness, persistence and source hypothesis assembled	Report assembly	<i>Investigation report</i>

STEPS IN WORKFLOW

18

TRAINING DURATION

86 h · 15 d

PROGRAMME CODE

SBQ-18

TRAINING FEE (INR)

2,00,000

Environmental & Utility Monitoring Metagenomics

(Facility environmental monitoring by sequencing · Environmental monitoring metagenomics)

WHY IT IS DONE

Culture-based environmental monitoring recovers only what grows on the chosen media, which is a small fraction of what is present. Sequencing characterises the facility's full microbial landscape, which is what makes trend analysis and source attribution possible.

WHAT IT ANSWERS

What microbial populations are present across this facility and its utilities, how do they change over time, and where do recurrent organisms persist?

WHEN IT IS DONE

In routine monitoring programmes alongside compendial methods, in facility qualification, in water system characterisation, and after a contamination event.

WHAT TRIGGERS IT

A routine monitoring cycle, facility or utility qualification, a contamination investigation, or a trend excursion.

HOW IT IS DONE

Sampling locations and controls are designed with the facility layout in mind, samples are processed with controls run identically because low-biomass environmental samples are dominated by reagent background. Communities are profiled taxonomically, background is subtracted, and results are analysed spatially and temporally against historical data.

WHAT IT PRODUCES

Community profiles by location with control-corrected results, organisms not recovered by culture, temporal trend analysis, spatial pattern mapping, recurrent organism identification, and monitoring programme recommendations.

WHAT IT DETECTS

Microbial community composition across sampled locations, organisms undetectable by the culture media in use, temporal trends and recurrent populations, water system biofilm signatures, and spatial patterns across the facility.

WHAT IT DOES NOT DETECT

It detects nucleic acid rather than viable organisms, so it does not replace compendial sterility or bioburden testing. Dead organisms and residual nucleic acid from cleaning are detected and can be misread as live contamination.

WHAT MAKES IT HARD

Environmental samples are low biomass and reagent background dominates them, so controls processed in parallel are the difference between a facility profile and a kit profile. Distinguishing residual nucleic acid from cleaned surfaces from live contamination is not possible by sequencing alone.

WHO DOES THIS JOB

Bioinformatician in quality microbiology, mid level, working with environmental monitoring programmes.

WHAT YOU CAN DO AFTERWARDS

Run facility environmental metagenomics with proper background control, analyse spatial and temporal patterns, and support monitoring programme design.

WHO HIRES FOR IT

Sterile and biologics manufacturers, contract manufacturing organisations, quality microbiology functions, and facility qualification consultancies.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
17	82 hours · 14 days	Prior metagenomic analysis experience	1,90,000

Environmental & Utility Monitoring Metagenomics

(Facility environmental monitoring by sequencing · Environmental monitoring metagenomics)

COMPLETE WORKFLOW — 17 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Monitoring programme scope confirmation	Locations, frequency and objectives of the monitoring programme confirmed	Programme review	<i>Programme scope</i>
02	Sampling design review	Sampling locations reviewed against facility layout and risk zones	Design review	<i>Sampling design</i>
03	Control design	Field, reagent and negative controls designed into every sampling round	Control design	<i>Control plan</i>
04	Sample and control processing	Samples and controls processed in parallel under identical conditions	Processing	<i>Prepared libraries</i>
05	Data receipt and quality assessment	Sequence integrity and quality verified per sample and control	md5sum, FastQC	<i>QC acceptance record</i>
06	Low-biomass handling verification	Processing confirmed appropriate for low-biomass environmental samples	Method verification	<i>Low-biomass record</i>
07	Taxonomic profiling	Community composition profiled per sampling location	Classification	<i>Community profiles</i>
08	Reagent background characterisation	Background taxa characterised from reagent and negative controls	Background analysis	<i>Background profile</i>
09	Background subtraction	Sample profiles corrected for reagent-derived background	Decontamination	<i>Corrected profiles</i>
10	Culture comparison	Sequencing results compared against parallel culture-based monitoring	Comparative analysis	<i>Culture comparison</i>
11	Non-culturable organism identification	Organisms undetected by the culture media in use identified	Comparative analysis	<i>Non-culturable findings</i>
12	Spatial pattern analysis	Community patterns analysed across facility zones and locations	Spatial analysis	<i>Spatial patterns</i>
13	Temporal trend analysis	Composition trends analysed against historical monitoring data	Trend analysis	<i>Temporal trends</i>
14	Recurrent organism identification	Organisms recurring across time and locations identified	Recurrence analysis	<i>Recurrent organisms</i>
15	Water system signature assessment	Biofilm and water system associated organisms assessed specifically	Targeted analysis	<i>Water system findings</i>
16	Viability limitation documentation	The nucleic acid versus viable organism distinction documented	Limitation template	<i>Viability caveat</i>
17	Programme recommendation	Monitoring programme and media recommendations derived from the findings	Recommendation framework	<i>Programme recommendations</i>

STEPS IN WORKFLOW

17

TRAINING DURATION

82 h · 14 d

PROGRAMME CODE

SBQ-19

TRAINING FEE (INR)

1,90,000

Process Comparability Genomic Analysis

(Manufacturing comparability genomic assessment · Comparability analysis)

WHY IT IS DONE

Every process change, site transfer and scale-up requires demonstration that the product is unchanged. Genomic attributes are increasingly part of that demonstration for vectors, cells and nucleic acid products, and comparability failures are among the costliest delays in manufacturing.

WHAT IT ANSWERS

Are the genomic quality attributes of product made by the new process equivalent to those made by the old one?

WHEN IT IS DONE

In process changes, site transfers, scale-up, and raw material or supplier changes affecting the product.

WHAT TRIGGERS IT

A planned process change, site transfer, scale-up, or a supplier change requiring comparability demonstration.

HOW IT IS DONE

Comparability attributes and acceptance criteria are pre-specified before any post-change material is analysed, historical batch variation is characterised to establish the comparison range, and post-change batches are analysed by identical methods. Differences are assessed statistically against historical variation, and any difference outside range is investigated for cause and consequence.

WHAT IT PRODUCES

Pre-specified comparability protocol with acceptance criteria, historical variation characterisation, attribute-by-attribute comparison with statistical assessment, difference investigation where applicable, and a comparability conclusion.

WHAT IT DETECTS

Differences in genomic quality attributes between pre-change and post-change material, whether differences fall within historical batch variation, attribute-by-attribute equivalence assessment, and trends indicating a systematic process shift.

WHAT IT DOES NOT DETECT

It does not establish clinical comparability, which is a broader assessment. Attributes not measured cannot be compared, so the comparability protocol determines what the exercise can conclude.

WHAT MAKES IT HARD

Historical variation must be characterised from enough batches to define a meaningful range, and most programmes have too few. Pre-specifying acceptance criteria before seeing post-change data is what makes the conclusion credible and is frequently reversed in practice.

WHO DOES THIS JOB

Bioinformatician in analytical development or quality, senior level, working with statistics and regulatory affairs.

WHAT YOU CAN DO AFTERWARDS

Design comparability protocols with pre-specified criteria, characterise historical variation, and assess attribute equivalence statistically.

WHO HIRES FOR IT

Biologics manufacturers, cell and gene therapy companies, contract manufacturers, and regulatory affairs functions.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
21	102 hours · 17 days	Prior analytical and statistical analysis experience	2,35,000

COMPLETE WORKFLOW — 21 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Change scope definition	The process change and its potential product impact defined	Change control record	<i>Change scope</i>
02	Comparability attribute selection	Genomic attributes relevant to the change selected and justified	Attribute selection	<i>Selected attributes</i>
03	Acceptance criteria pre-specification	Acceptance criteria specified in writing before post-change analysis	Protocol approval	<i>Pre-specified criteria</i>
04	Historical batch identification	Pre-change batches identified for characterising historical variation	Batch selection	<i>Historical batch set</i>
05	Method consistency verification	Analytical methods confirmed identical across pre and post-change material	Method verification	<i>Method consistency record</i>
06	Historical variation characterisation	Batch-to-batch variation characterised from pre-change material	Statistical analysis	<i>Historical variation</i>
07	Variation adequacy assessment	Whether enough batches exist to define a meaningful range assessed	Adequacy assessment	<i>Adequacy determination</i>
08	Data receipt and quality assessment	Data quality verified for all batches under identical criteria	Quality assessment	<i>QC acceptance record</i>
09	Post-change batch analysis	Post-change batches analysed by the identical analytical method	Analysis execution	<i>Post-change results</i>
10	Attribute-by-attribute comparison	Each attribute compared between pre and post-change material	Comparative analysis	<i>Attribute comparisons</i>
11	Statistical equivalence assessment	Differences assessed statistically against historical variation	Statistical testing	<i>Equivalence assessment</i>
12	Acceptance criteria evaluation	Results evaluated against the pre-specified acceptance criteria	Criteria evaluation	<i>Criteria evaluation</i>
13	Out-of-range difference investigation	Differences outside historical range investigated for cause	Investigation procedure	<i>Difference investigation</i>
14	Consequence assessment	Any confirmed difference assessed for product quality consequence	Assessment framework	<i>Consequence assessment</i>
15	Trend assessment	Post-change batches assessed for systematic trend rather than scatter	Trend analysis	<i>Trend assessment</i>
16	Additional testing determination	Whether further characterisation is required determined	Determination criteria	<i>Additional testing decision</i>
17	Comparability conclusion	A comparability conclusion drawn against the pre-specified criteria	Conclusion framework	<i>Comparability conclusion</i>
18	Risk assessment	Residual risk from any unresolved difference assessed	Risk assessment	<i>Risk determination</i>
19	Regulatory strategy support	Findings assembled to support the regulatory strategy for the change	Support assembly	<i>Regulatory support</i>
20	Comparability report assembly	Criteria, variation, comparisons and conclusion assembled	Report assembly	<i>Comparability report</i>
21	Change control closure support	Evidence assembled to support closure of the change control record	Documentation	<i>Closure support</i>

STEPS IN WORKFLOW

21

TRAINING DURATION

102 h · 17 d

PROGRAMME CODE

SBQ-20

TRAINING FEE (INR)

2,35,000

Stability Study Genomic Analysis

(Product stability genomic assessment · Stability genomic analysis)

WHY IT IS DONE

Shelf life is set by the attribute that degrades first, and for nucleic acid and vector products that is frequently a genomic attribute rather than a protein one. Stability data determine storage conditions, shelf life and shipping constraints.

WHAT IT ANSWERS

Do the genomic quality attributes of this product change over storage, at what rate, and what shelf life do they support?

WHEN IT IS DONE

In formal stability programmes at defined timepoints, in accelerated and stress studies, and in shipping and handling qualification.

WHAT TRIGGERS IT

A stability programme timepoint, an accelerated study, a shipping qualification, or an excursion investigation.

HOW IT IS DONE

Stability timepoints and conditions are executed with identical analytical methods throughout, attributes are measured at each timepoint, and trends are modelled rather than assessed pairwise. Degradation rates are estimated, shelf life is projected from the trend, and accelerated and real-time data are compared for pathway consistency.

WHAT IT PRODUCES

Attribute measurements across timepoints and conditions, trend models with degradation rate estimates, shelf life projection with confidence, accelerated against real-time comparison, and storage condition recommendations.

WHAT IT DETECTS

Change in sequence integrity over time, fragmentation and degradation trends, impurity species increase, attribute trajectories against time and condition, and the timepoint at which any attribute approaches specification.

WHAT IT DOES NOT DETECT

It does not measure potency or biological activity directly, and genomic degradation does not always correlate with potency loss. Accelerated conditions may produce degradation pathways that differ from real-time storage.

WHAT MAKES IT HARD

Method drift across a multi-year study is indistinguishable from product degradation unless controls are analysed alongside every timepoint. Trend modelling rather than pairwise comparison is what supports a shelf life claim, and pairwise testing is what most studies actually do.

WHO DOES THIS JOB

Bioinformatician in analytical development, mid level, working with stability programme management.

WHAT YOU CAN DO AFTERWARDS

Analyse stability data as trends, control for method drift, estimate degradation rates, and support shelf life determination.

WHO HIRES FOR IT

Biologics, vector and mRNA manufacturers, analytical development functions, and contract testing organisations.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
16	78 hours · 13 days	Prior analytical development experience	1,80,000

COMPLETE WORKFLOW — 16 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Stability protocol confirmation	Timepoints, conditions and attributes confirmed against the approved protocol	Protocol review	Stability protocol
02	Method consistency planning	Analytical method consistency across the study duration planned	Method planning	Consistency plan
03	Control material designation	Control material designated for analysis alongside every timepoint	Control designation	Study controls
04	Data receipt and quality assessment	Data quality verified per timepoint against identical criteria	Quality assessment	QC acceptance record
05	Control result review	Control results reviewed at each timepoint to detect method drift	Control analysis	Method drift assessment
06	Attribute measurement	Each stability attribute measured at every timepoint and condition	Analysis execution	Timepoint measurements
07	Method drift correction	Measurements corrected or qualified where control drift is detected	Correction procedure	Corrected measurements
08	Trend modelling	Attributes modelled as trends across time rather than compared pairwise	Trend modelling	Trend models
09	Degradation rate estimation	Degradation rates estimated per attribute and condition	Rate estimation	Degradation rates
10	Condition comparison	Degradation compared across storage conditions	Comparative analysis	Condition comparison
11	Accelerated to real-time comparison	Accelerated and real-time degradation pathways compared for consistency	Pathway comparison	Pathway consistency
12	Specification approach projection	The timepoint at which each attribute approaches specification projected	Projection modelling	Approach projection
13	Shelf life derivation	Shelf life derived from the limiting attribute with confidence stated	Derivation framework	Shelf life projection
14	Storage condition recommendation	Storage and handling conditions recommended from the findings	Recommendation framework	Storage recommendations
15	Excursion impact assessment	Impact of temperature or handling excursions assessed where relevant	Impact assessment	Excursion assessment
16	Stability report assembly	Measurements, trends, rates and shelf life projection assembled	Report assembly	Stability report

STEPS IN WORKFLOW	TRAINING DURATION	PROGRAMME CODE	TRAINING FEE (INR)
16	78 h · 13 d	SBQ-21	1,80,000

Analytical Method Validation for Sequencing Assays

(Sequencing assay validation for GMP use · Method validation)

WHY IT IS DONE

A sequencing method used in a regulated release decision must be validated to the same standard as any other analytical method, and the validation parameters that apply to chemistry do not translate straightforwardly to sequencing. Getting this right is what allows a method to be used at all.

WHAT IT ANSWERS

Does this sequencing method meet the validation requirements for its intended use, across every applicable parameter?

WHEN IT IS DONE

Before any sequencing method is used for release, comparability or stability decisions, and when a validated method is changed.

WHAT TRIGGERS IT

A method progressing to regulated use, a method change requiring revalidation, or a regulatory question on method suitability.

HOW IT IS DONE

Intended use and required performance are defined first because they determine every study, reference materials are selected, and each validation parameter is executed as a designed study rather than an observation. Bioinformatic pipeline contribution is validated separately from wet-lab performance, robustness is tested by deliberate variation, and the validated operating range is defined explicitly.

WHAT IT PRODUCES

Intended use and requirement definition, accuracy, precision, specificity, detection limit, linearity and robustness study results, separate pipeline validation evidence, defined operating range, and a validation report with acceptance against criteria.

WHAT IT DETECTS

Method accuracy against reference material, precision within and between runs and operators, specificity and interference, limits of detection and quantitation, linearity and range, robustness to deliberate variation, and the boundaries of the validated method.

WHAT IT DOES NOT DETECT

Validation establishes that the method performs as characterised, not that the attribute it measures is the right one. It also cannot anticipate every sample condition encountered in routine use.

WHAT MAKES IT HARD

The bioinformatic pipeline is part of the method and must be validated as such, including version control and change impact, which is where most sequencing method validations are weakest. Detection limit in sequencing depends on depth and background and must be established per matrix rather than claimed generally.

WHO DOES THIS JOB

Bioinformatician or analytical scientist in method validation, senior level, working with quality assurance and regulatory affairs.

WHAT YOU CAN DO AFTERWARDS

Validate sequencing methods to regulated standards, validate the analytical pipeline as part of the method, and define operating ranges defensibly.

WHO HIRES FOR IT

Biologics and advanced therapy manufacturers, contract testing laboratories, analytical development functions, and regulatory consultancies.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
24	118 hours · 20 days	Prior assay development and statistical experience	2,70,000

Analytical Method Validation for Sequencing Assays

(Sequencing assay validation for GMP use · Method validation)

COMPLETE WORKFLOW — 24 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Intended use definition	The intended use of the method defined since it determines every requirement	Use specification	<i>Intended use statement</i>
02	Regulatory expectation review	Applicable validation expectations for the method type reviewed	Guidance review	<i>Expectation record</i>
03	Performance requirement derivation	Required performance derived from intended use and specification limits	Requirement derivation	<i>Performance requirements</i>
04	Method description documentation	The complete method including wet-lab and computational steps documented	Method documentation	<i>Method description</i>
05	Reference material selection	Reference and control materials selected and characterised	Material selection	<i>Reference materials</i>
06	Validation protocol authorisation	The validation protocol and acceptance criteria approved before execution	Protocol approval	<i>Authorised protocol</i>
07	Specificity study execution	Specificity established including interference and cross-reactivity	Specificity study	<i>Specificity results</i>
08	Accuracy study execution	Accuracy established against characterised reference material	Accuracy study	<i>Accuracy results</i>
09	Repeatability study execution	Within-run precision established across replicates	Repeatability study	<i>Repeatability results</i>
10	Intermediate precision study execution	Between-run, operator and instrument precision established	Precision study	<i>Intermediate precision</i>
11	Detection limit determination	Limit of detection established per matrix rather than claimed generally	Detection limit study	<i>Detection limits</i>
12	Quantitation limit determination	Limit of quantitation established where the method is quantitative	Quantitation study	<i>Quantitation limits</i>
13	Linearity and range study	Linearity and the working range established	Linearity study	<i>Linearity results</i>
14	Robustness study execution	Robustness tested by deliberate variation of critical parameters	Robustness study	<i>Robustness results</i>
15	Computational pipeline validation	The analytical pipeline validated independently of wet-lab performance	Pipeline validation	<i>Pipeline validation record</i>
16	Pipeline version control establishment	Version control and change control established for the pipeline	Version control	<i>Version control record</i>
17	Reference data reproducibility testing	Reference datasets re-analysed to confirm pipeline reproducibility	Reproducibility testing	<i>Reproducibility results</i>
18	Sample stability assessment	Sample and extract stability under handling conditions assessed	Stability study	<i>Sample stability</i>
19	Carryover and contamination assessment	Run-to-run carryover and cross-contamination assessed	Carryover study	<i>Carryover results</i>
20	Operating range definition	The validated operating range defined explicitly with boundaries stated	Range definition	<i>Operating range</i>
21	Acceptance criteria evaluation	All results evaluated against pre-specified acceptance criteria	Criteria evaluation	<i>Acceptance evaluation</i>
22	System suitability criteria derivation	Routine system suitability criteria derived from validation data	Derivation framework	<i>System suitability criteria</i>
23	Revalidation trigger definition	Conditions requiring revalidation defined and documented	Trigger definition	<i>Revalidation triggers</i>
24	Validation report assembly	All studies, results and conclusions assembled into a validation report	Report assembly	<i>Validation report</i>

STEPS IN WORKFLOW

24

TRAINING DURATION

118 h · 20 d

PROGRAMME CODE

SBQ-22

TRAINING FEE (INR)

2,70,000

Data Integrity & Computerised System Compliance Analysis

(Data integrity assessment for computational systems · Data integrity compliance)

WHY IT IS DONE

Regulatory findings against manufacturers are dominated by data integrity rather than science, and computational systems handling sequencing data are increasingly within inspection scope. A correct result that cannot be shown to be attributable, legible and contemporaneous is a finding regardless of its correctness.

WHAT IT ANSWERS

Do the computational systems producing this data meet data integrity expectations, and can every reported result be traced and reproduced?

WHEN IT IS DONE

In inspection readiness preparation, in qualifying a new computational system, after a data integrity observation, and in periodic system review.

WHAT TRIGGERS IT

Inspection readiness, new system qualification, a regulatory observation, or a scheduled periodic review.

HOW IT IS DONE

System scope and the records within it are defined, audit trail coverage and review practice are assessed, and traceability is tested by selecting reported results and tracing them back to raw data. Access control, change control and backup arrangements are reviewed, reproducibility is tested by re-execution, and gaps are documented with remediation.

WHAT IT PRODUCES

System scope and record inventory, audit trail coverage assessment, traceability test results, access and change control review, backup and retention assessment, reproducibility test results, and a gap and remediation register.

WHAT IT DETECTS

Gaps in audit trail coverage, breaks in traceability from raw data to reported result, uncontrolled analytical changes, access control weaknesses, backup and retention gaps, and results that cannot be reproduced from retained data.

WHAT IT DOES NOT DETECT

It assesses control and documentation rather than scientific correctness. A fully compliant system can produce a scientifically wrong answer, and that is a separate review.

WHAT MAKES IT HARD

Retrospective assessment finds gaps that cannot be closed retrospectively, because contemporaneous records cannot be recreated. Analytical scripts and manual data handling steps outside validated systems are where traceability usually breaks and are the hardest to bring under control.

WHO DOES THIS JOB

Bioinformatician or computational quality specialist, senior level, working with quality assurance and information technology.

WHAT YOU CAN DO AFTERWARDS

Assess computational systems against data integrity expectations, test traceability and reproducibility, and construct remediation plans.

WHO HIRES FOR IT

Biologics and pharmaceutical manufacturers, contract organisations, quality assurance functions, and compliance consultancies.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
20	98 hours · 17 days	Prior pipeline and quality systems experience	2,25,000

Data Integrity & Computerised System Compliance Analysis

(Data integrity assessment for computational systems · Data integrity compliance)

COMPLETE WORKFLOW — 20 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	System scope definition	Computational systems and records within scope identified	Scope definition	<i>System scope</i>
02	Record inventory	Electronic records produced by each system inventoried and classified	Inventory compilation	<i>Record inventory</i>
03	Regulatory expectation mapping	Applicable data integrity expectations mapped to each record type	Expectation mapping	<i>Expectation map</i>
04	Audit trail coverage assessment	Audit trail coverage assessed against the records requiring it	Coverage assessment	<i>Audit trail assessment</i>
05	Audit trail review practice assessment	Whether audit trails are reviewed and by whom assessed	Practice review	<i>Review practice record</i>
06	Access control review	User access, privileges and segregation of duties reviewed	Access review	<i>Access control assessment</i>
07	Change control assessment	Change control over analytical software and configuration assessed	Change control review	<i>Change control assessment</i>
08	Version control assessment	Version control over pipelines, scripts and references assessed	Version control review	<i>Version control assessment</i>
09	Traceability test sampling	Reported results sampled for end-to-end traceability testing	Sampling procedure	<i>Traceability sample</i>
10	Raw data traceability testing	Each sampled result traced back to its raw data through every step	Traceability testing	<i>Traceability results</i>
11	Intermediate step documentation assessment	Manual steps and scripts outside validated systems identified	Gap analysis	<i>Manual step inventory</i>
12	Reproducibility testing	Sampled results re-executed to confirm reproducibility	Re-execution	<i>Reproducibility results</i>
13	Discrepancy investigation	Any failure to reproduce investigated and explained	Investigation procedure	<i>Discrepancy record</i>
14	Backup and retention assessment	Backup, archival and retention arrangements assessed against requirements	Retention review	<i>Retention assessment</i>
15	Restoration testing	Restoration from backup tested and verified	Restoration test	<i>Restoration results</i>
16	Data lifecycle assessment	The full data lifecycle from generation to disposal assessed	Lifecycle review	<i>Lifecycle assessment</i>
17	Gap identification and classification	Identified gaps classified by severity and regulatory exposure	Gap classification	<i>Gap register</i>
18	Remediation planning	Remediation actions planned with owners and timelines	Remediation planning	<i>Remediation plan</i>
19	Retrospective limitation documentation	Gaps that cannot be closed retrospectively documented explicitly	Limitation documentation	<i>Retrospective limitations</i>
20	Assessment report assembly	Scope, findings, gaps and remediation assembled	Report assembly	<i>Data integrity report</i>

STEPS IN WORKFLOW

20

TRAINING DURATION

98 h · 17 d

PROGRAMME CODE

SBQ-23

TRAINING FEE (INR)

2,25,000

Deviation & Out-of-Specification Investigation Analysis

(Genomic investigation of manufacturing deviations · Deviation investigation)

WHY IT IS DONE

When a batch fails specification the investigation determines whether it is released, reworked or discarded, and whether the process is fit to continue. Genomic evidence frequently distinguishes a laboratory error from a genuine product problem, which are handled entirely differently.

WHAT IT ANSWERS

What caused this out-of-specification or deviation result, was it laboratory error or a genuine product attribute, and what is the impact on this and other batches?

WHEN IT IS DONE

Immediately following an out-of-specification result, in deviation investigation, and in assessing impact across other batches.

WHAT TRIGGERS IT

An out-of-specification result, a manufacturing deviation, an unexpected trend, or a complaint investigation.

HOW IT IS DONE

The original result is assessed for analytical validity before any retesting, sample identity is verified to exclude mix-up, and the analysis is repeated under controlled conditions. Contamination and error signatures are sought specifically, candidate root causes are tested against the evidence, and impact on other batches is assessed by retrospective analysis.

WHAT IT PRODUCES

Original result validity assessment, sample identity verification, controlled retest results, error signature analysis, root cause evidence assessment, batch impact analysis, and an investigation conclusion supporting disposition.

WHAT IT DETECTS

Whether the original result is reproducible, laboratory error signatures such as sample mix-up or contamination, genuine product attribute deviation, the extent of impact across batches, and evidence supporting or refuting candidate root causes.

WHAT IT DOES NOT DETECT

It cannot establish root cause on its own, since the analytical evidence must be combined with process and laboratory records. An unreproducible result does not prove laboratory error, only that the original observation could not be confirmed.

WHAT MAKES IT HARD

Retesting without first establishing analytical validity of the original result is the classic investigation failure, because it converts an investigation into a search for a passing result. Sample identity verification is what most often resolves an investigation and is most often done last.

WHO DOES THIS JOB

Bioinformatician in quality investigation, mid to senior level, working with quality assurance and manufacturing.

WHAT YOU CAN DO AFTERWARDS

Investigate out-of-specification results systematically, distinguish laboratory error from product deviation, and assess batch impact.

WHO HIRES FOR IT

Biologics and advanced therapy manufacturers, contract manufacturing organisations, quality assurance functions, and quality control laboratories.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
19	92 hours · 16 days	Prior analytical and investigative experience	2,10,000

Deviation & Out-of-Specification Investigation Analysis

(Genomic investigation of manufacturing deviations · Deviation investigation)

COMPLETE WORKFLOW — 19 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Event capture and scope definition	The deviation or out-of-specification result captured with full context	Event record	<i>Event record</i>
02	Immediate containment assessment	Whether affected material requires immediate quarantine assessed	Containment procedure	<i>Containment decision</i>
03	Original result validity assessment	The original result assessed for analytical validity before any retesting	Validity assessment	<i>Validity determination</i>
04	Analytical record review	Instrument, run and analytical records reviewed for the original result	Record review	<i>Record review findings</i>
05	Sample identity verification	Sample identity verified to exclude mix-up as the cause	Identity analysis	<i>Identity verification</i>
06	Sample integrity assessment	Sample handling and integrity assessed as a candidate cause	Integrity assessment	<i>Integrity findings</i>
07	Contamination signature analysis	The result assessed for signatures characteristic of contamination	Signature analysis	<i>Contamination assessment</i>
08	Laboratory error hypothesis testing	Laboratory error hypotheses tested against the available evidence	Hypothesis testing	<i>Error hypothesis results</i>
09	Controlled retest design	Retesting designed with controls and criteria specified in advance	Retest design	<i>Retest protocol</i>
10	Controlled retest execution	Retesting executed under the specified protocol	Retest execution	<i>Retest results</i>
11	Retest interpretation	Retest results interpreted against the original without bias toward passing	Interpretation framework	<i>Retest interpretation</i>
12	Product attribute hypothesis testing	Genuine product deviation hypotheses tested against the evidence	Hypothesis testing	<i>Product hypothesis results</i>
13	Process record correlation	Findings correlated against manufacturing and process records	Correlation analysis	<i>Process correlation</i>
14	Root cause evidence assessment	Candidate root causes assessed against all assembled evidence	Assessment framework	<i>Root cause assessment</i>
15	Batch impact analysis	Impact on other batches assessed by retrospective analysis	Retrospective analysis	<i>Batch impact</i>
16	Trend review	Historical data reviewed for prior occurrences of the same signal	Trend analysis	<i>Trend findings</i>
17	Corrective action support	Evidence assembled to support corrective and preventive action	Support assembly	<i>Corrective action support</i>
18	Disposition recommendation	A disposition recommendation made for the affected material	Recommendation framework	<i>Disposition recommendation</i>
19	Investigation report assembly	Validity, testing, root cause and impact assembled into the investigation record	Report assembly	<i>Investigation report</i>

STEPS IN WORKFLOW

19

TRAINING DURATION

92 h · 16 d

PROGRAMME CODE

SBQ-24

TRAINING FEE (INR)

2,10,000

Rapid Sterility Assurance by Sequencing Methods

(Rapid microbiological method implementation · Rapid sterility methods)

WHY IT IS DONE

Compendial sterility testing takes fourteen days, which is longer than the shelf life of a cell therapy product. Sequencing-based rapid methods provide an answer within the release window, and demonstrating their equivalence is what allows them to be used.

WHAT IT ANSWERS

Does this product show microbial sequence within the release window, and is the rapid method demonstrably equivalent to the compendial one?

WHEN IT IS DONE

In release testing of short-shelf-life products, in method implementation and equivalence demonstration, and in supporting regulatory acceptance of a rapid method.

WHAT TRIGGERS IT

A short-shelf-life product requiring release, a rapid method implementation project, or a regulatory submission for method acceptance.

HOW IT IS DONE

The method is challenged against the required organism panel at low inoculum with detection limits established per organism, and equivalence against the compendial method is demonstrated across the panel. Routine samples are analysed with controls processed identically, detections are confirmed and identified, and results are reported within the release window.

WHAT IT PRODUCES

Per-organism detection limits from challenge studies, compendial equivalence demonstration, routine sample results with control comparison, detection identification and confirmation, and a release recommendation within the window.

WHAT IT DETECTS

Microbial sequence at the sensitivity the method is validated to, species identification of detections, and the achieved detection limit for the panel of challenge organisms.

WHAT IT DOES NOT DETECT

It detects nucleic acid rather than viable organisms, which is the central obstacle to accepting it as a sterility method, and dead organisms from upstream processing are detected. Very low-level contamination near a single organism per container remains at the edge of what any method resolves.

WHAT MAKES IT HARD

Demonstrating equivalence to a fourteen-day culture method requires challenge across a defined organism panel at very low inoculum, and sensitivity differs substantially between organisms. The viable-versus-detectable distinction is the argument regulators press hardest on.

WHO DOES THIS JOB

Bioinformatician in quality microbiology or rapid methods, senior level, working with microbiology and regulatory affairs.

WHAT YOU CAN DO AFTERWARDS

Implement and validate rapid sequencing-based microbiological methods, demonstrate compendial equivalence, and operate them within release windows.

WHO HIRES FOR IT

Cell and gene therapy manufacturers, sterile product manufacturers, rapid method technology providers, and quality microbiology functions.

COMMERCIALS

STEPS IN WORKFLOW	TRAINING DURATION	PREREQUISITE	TRAINING FEE (INR)
18	88 hours · 15 days	Prior microbial detection analysis experience	2,05,000

Rapid Sterility Assurance by Sequencing Methods

(Rapid microbiological method implementation · Rapid sterility methods)

COMPLETE WORKFLOW — 18 ORDERED STEPS

#	STEP	WHAT HAPPENS	TOOL / METHOD	OUTPUT
01	Method scope and intended use definition	The products and release decisions the method will support defined	Use specification	<i>Intended use statement</i>
02	Compendial comparison basis establishment	The compendial method against which equivalence will be shown identified	Comparison basis	<i>Comparison basis</i>
03	Challenge organism panel definition	The organism panel for challenge studies defined against requirements	Panel definition	<i>Challenge panel</i>
04	Low inoculum challenge preparation	Challenge organisms prepared at the low inoculum levels required	Inoculum preparation	<i>Challenge inocula</i>
05	Matrix-specific challenge execution	Challenges performed in the actual product matrix	Challenge execution	<i>Challenge results</i>
06	Per-organism detection limit determination	Detection limits established for every organism in the panel	Limit determination	<i>Per-organism limits</i>
07	Compendial parallel testing	Compendial testing performed in parallel on the same challenged material	Parallel testing	<i>Compendial results</i>
08	Equivalence assessment	Rapid method performance assessed against compendial across the panel	Equivalence analysis	<i>Equivalence assessment</i>
09	Inhibition and interference assessment	Product matrix inhibition of detection assessed	Interference study	<i>Interference assessment</i>
10	Control design for routine use	Positive, negative and process controls designed for routine operation	Control design	<i>Routine control plan</i>
11	Routine sample analysis	Product samples analysed with controls processed in parallel	Analysis execution	<i>Sample results</i>
12	Detection assessment against controls	Detections assessed against parallel controls before interpretation	Control comparison	<i>Detection assessment</i>
13	Organism identification	Detected organisms identified to species where possible	Identification analysis	<i>Identification results</i>
14	Confirmation analysis	Detections confirmed by an orthogonal or targeted approach	Confirmation analysis	<i>Confirmation results</i>
15	Viability limitation documentation	The viable versus detectable distinction documented explicitly	Limitation template	<i>Viability caveat</i>
16	Turnaround verification	Result delivery verified against the release window requirement	Timing audit	<i>Turnaround record</i>
17	Release recommendation	A release recommendation made within the available window	Recommendation framework	<i>Release recommendation</i>
18	Method report assembly	Challenge results, equivalence, detections and recommendation assembled	Report assembly	<i>Rapid method report</i>

STEPS IN WORKFLOW

18

TRAINING DURATION

88 h · 15 d

PROGRAMME CODE

SBQ-25

TRAINING FEE (INR)

2,05,000

Build and staff a *GMP genomic quality function*

Each of the twenty-five programmes in this catalogue trains one complete type of NGS data analysis as a regulated manufacturing quality function performs it — not a tool, not a step, but the whole analysis from the data that arrives to the report that is signed. Trainees run the pipeline on live clinical data, break it, diagnose the failure, and produce the deliverable report set a laboratory would place in a case file and defend to an assessor. Programmes are named the way laboratories name their own work, so a trainee returns to the bench able to run the analysis on the requisition rather than a generic version of it. Individual programmes may be taken alone, combined into a track, or delivered in-house at the client laboratory against that laboratory's own assays and gene lists.

CATALOGUE

SBQ · Sector 05 of 16

DURATION

2328 contact hours · 397 working days

SCOPE

25 analysis types · 481 steps

COMPLETE CATALOGUE FEE

INR 34,84,000 (list 53,60,000)

15%

ANY THREE
PROGRAMMES

22%

CORE
TRACK

28%

PROFESSIONAL
TRACK

35%

COMPLETE
CATALOGUE